

Certificate of Conformity Finntip™ Filter Sterile

Product Code: 94052310
Product Name: Finntip Filter 100 µL Ext, sterile
Lot Number: 19174C0
Exp: 2024-06

We hereby confirm that these high quality Finntip™ Pipette tips are manufactured under controlled environment and are certified by validation of the manufacturing process to be free of the following contaminants:

| | |
|-----------|-------------------------------|
| Human DNA | < 1.9×10^{-11} g/tip |
| DNase | < 9.4×10^{-4} U/tip |
| RNase | < 6.5×10^{-8} U/tip |
| Endotoxin | < 0.05 EU/ml |

The sterility (SAL 10^{-6}) is ensured by irradiating this product according to ISO 11137 guidelines.

In Joensuu, 2nd December 2019



Heidi Salokangas, Quality Manager
Laboratory Plastics Essentials
Thermo Fisher Scientific Oy