

VWR® STORAGE CARTS

100% Total Replacement
Lifetime Guarantee

Available in 7 Colors

Optional Railtop

Chemically Resistant
Polyethylene
Construction




VWR® STORAGE CARTS



Dispensing and Regular Side Storage Shelves



Bulk Storage Area with Magnetic Door




Lock Drawers & Bulk Storage with 1 Key



Full-Extension Drawers with Handle



Protective Film Top




Fully Enclosed 3" Casters


Select the Railed Top for Additional Security

CHOOSE FROM 7 UNIQUE COLORS

WHITE




BLACK



METALLIC SILVER




ALMOND BEIGE



GUM DROP PURPLE



BAHAMAS SEA TEAL



GLOBAL BLUE



Need another color? Check out the back page to personalize even more!

Mobilize. Organize. Personalize.

5 Sizes



	Mini	Small	Medium	Large	Jumbo
Description	PL22 Mini Storage Cart	PL27 Small Storage Cart	PL31 Medium Storage Cart	PL36 Large Storage Cart	PL41 Jumbo Storage Cart
External Dimensions (W x H x D)	22 ⁷ / ₈ x 22 x 18 ¹ / ₄ "	22 ⁷ / ₈ x 27 x 18 ¹ / ₄ "	22 ⁷ / ₈ x 30 ⁵ / ₈ x 18 ¹ / ₄ "	22 ⁷ / ₈ x 35 ⁵ / ₈ x 18 ¹ / ₄ "	22 ⁷ / ₈ x 40 ¹ / ₂ x 18 ¹ / ₄ "
Number of Drawers	1 Drawer	2 Drawers	2 Drawers	3 Drawers	4 Drawers
Bulk Storage Internal Dimensions (W x H x D)	12 ⁷ / ₈ x 11 ³ / ₈ x 15 ¹ / ₈ "	12 ⁷ / ₈ x 11 ³ / ₈ x 15 ¹ / ₈ "	12 ⁷ / ₈ x 15 ¹ / ₈ x 15 ¹ / ₈ "	12 ⁷ / ₈ x 15 ¹ / ₈ x 15 ¹ / ₈ "	12 ⁷ / ₈ x 15 ¹ / ₈ x 15 ¹ / ₈ "
VWR® Storage Carts with Flat Top					
White	10809-048	10809-050	10809-052	10809-054	10809-056
Global Blue	75857-056	75857-070	75857-084	75857-098	75857-112
Metallic Silver	75857-058	75857-072	75857-086	75857-100	75857-114
Almond Beige	75857-060	75857-074	75857-088	75857-102	75857-116
Gum Drop Purple	75857-064	75857-078	75857-092	75857-106	75857-120
Bahamas Sea Teal	75857-066	75857-080	75857-094	75857-108	75857-122
Black	75857-062	75857-076	75857-090	75857-104	75857-118
VWR® Storage Carts with Railed Top					
White	75857-126	75857-128	75857-130	75857-132	75857-134
Global Blue	75870-344	75870-342	75870-364	75870-378	75870-392
Metallic Silver	75870-346	75870-352	75870-366	75870-380	75870-394
Almond Beige	75870-348	75870-354	75870-368	75870-382	75870-396
Gum Drop Purple	75870-336	75870-358	75870-372	75870-386	75870-400
Bahamas Sea Teal	75870-338	75870-360	75870-374	75870-388	75870-402
Black	75870-350	75870-356	75870-370	75870-384	75870-398

All Drawers are the Same Size Internal Drawer Dimensions: 11 W x 3³/₈ H x 15¹/₈ " D




MOBILIZE. ORGANIZE. PERSONALIZE.



Customize your VWR® STORAGE CART

Add your organization's logo and/or custom drawer color



Just provide the logo artwork (vector format) or color sample

	Mini	Small	Medium	Large	Jumbo
Description	PL22 Mini Storage Cart	PL27 Small Storage Cart	PL31 Medium Storage Cart	PL36 Large Storage Cart	PL41 Jumbo Storage Cart
VWR® Storage Carts with Flat Top					
Customized	75857-068	75857-082	75857-096	75857-110	75857-124
VWR® Storage Carts with Railed Top					
Customized	75870-340	75870-362	75870-376	75870-390	75870-404



1.800.932.5000 | vwr.com

8044_RB_JS_2017

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.