



Item # _____

Job _____

Built for Labs and Cleanrooms using the most stringent cleaning agents.

316 Stainless Steel Table Top

WORK CENTERS

- **316 Table Tops** can withstand corrosive attacks by sodium, calcium brines, hypochlorite solutions, phosphoric acids, sulfate and sulfuric acid.
- **Built to Last:** The Space Saver 24" (610mm) Deep Stationary Worktables' sturdy construction provides a robust 800 lb. (363kg) load rating. Mobile tables are rated to 600 lbs.(272kg). Available preconfigured units consist of T-316 stainless steel components that stand up to the harshest applications and environments.
- **Proven Stability:** Space Saver worktables feature 1" (25mm) diameter posts and the tapered sleeves and corners of patented Super Erecta Shelf® construction. These features assure a rigid, sway-free table. The work surface is made of seamless, heavy-duty 16-gauge stainless steel that virtually eliminates any "drumming."
- **Easy to Assemble:** Units are shipped knocked down to save on freight, and can be put together in minutes without the use of special tools. Bottom shelf or frame is adjustable in 1" (25mm) increments. Table top fits over split sleeves at the top of the posts.
- **Stationary or Mobile:** Preconfigured units are offered both in stationary with a leveling foot and mobile utilizing a 5" (127mm) polymer caster. These innovative stem casters offer corrosion resistance and enhanced durability.
- **Popular Sizes:** Space Saver Worktables are offered in 24" (610mm) depth and 36" (914mm), 48" (1219mm) and 60" (1524mm) lengths. All models are 34" (860mm) high.
- **Customizable:** Preconfigured units are available in our T-316 stainless steel finish. Worktables can be ordered by component using Super Erecta® posts, frames and shelves in any finish for a solution that fits any price point or application.



*Stationary Worktable
with 3-sided
Double Snake Frame*



*Mobile Worktable with
Stainless Steel
Solid Undershaft*



InterMetro Industries Corporation
North Washington Street
Wilkes-Barre, PA 18705
www.metro.com



Space Saver 24" (610mm) Deep Worktables

17.42

Job _____

SPACE SAVER



24" (610mm) DEEP WORKTABLES

316 Stainless Steel Table Top

Specifications:

- **Top:** Seamless 16 gauge T-316 stainless steel with 11 gauge T-316 stainless steel channel support.
- **Bottom Shelf:** 18 gauge T-304 stainless steel with black, epoxy-coated corners. 36" length worktables utilize 21x30" (533x762mm) undershelves, 48" length worktables utilize 21x42" (533x1067mm) undershelves, and 60" length worktables utilize 21x48" (533x1219mm) undershelves. Can be used with any Super Erecta® or Super Adjustable Super Erecta® wire or solid shelving.
- **Bottom Frame:** 3-sided double snake frame offered in 304 stainless steel. 36" length worktables utilize 21x30" (533x762mm) frames, 48" length worktables utilize 21x42" (533x1067mm) frames, and 60" length worktables utilize 21x48" (533x1219mm) frames.
- **Posts:** Stainless Steel posts have an outside diameter of 1" (25mm), with rolled grooves in 1" (25mm) increments along entire height.
- **Casters:** Mobile units feature 2 swivel/2 brake 5" (127mm) polymer caster with polyurethane tread. Can also be used with any standard Metro Stem Caster. (Note: changing casters may affect load rating. See Specification Sheet 11.20 for more details)

Dimensions

Space Saver Worktables

24" (610mm) Wide

Stationary

Model No.	Length		Bottom Type	Approx. Pkd. Wt.	
	(in.)	(mm)		(lbs.)	(kg)
WTS2436FS-316T	36	914	Solid Stainless Steel	59.94	27.18
WTS2448FS-316T	48	1219	Solid Stainless Steel	77.60	35.20
WTS2460FS-316T	60	1524	Solid Stainless Steel	89.44	40.57
WTS2436US-316T	36	914	3-Sided Frame	51.44	23.33
WTS2448US-316T	48	1219	3-Sided Frame	63.04	28.59
WTS2460US-316T	60	1524	3-Sided Frame	73.74	33.45

Load Rating: Stationary Units are rated to 800 lbs. (363kg).

Mobile

Model No.	Length		Bottom Type	Approx. Pkd. Wt.	
	(in.)	(mm)		(lbs.)	(kg)
MWTS2436FS-316T	36	914	Solid Stainless Steel	65.02	29.49
MWTS2448FS-316T	48	1219	Solid Stainless Steel	82.68	37.50
MWTS2460FS-316T	60	1524	Solid Stainless Steel	94.52	42.87
MWTS2436US-316T	36	914	3-Sided Frame	56.52	25.64
MWTS2448US-316T	48	1219	3-Sided Frame	68.12	30.90
MWTS2460US-316T	60	1524	3-Sided Frame	78.82	35.75

Load Rating: Mobile Units are rated to 600 lbs. (272kg), evenly distributed. Mobile tables are intended for move-to-clean applications only and not as a transport device.All Metro Catalog Sheets are available on our website: www.metro.com

InterMetro Industries Corporation

North Washington Street, Wilkes-Barre, PA 18705

Phone: 570-825-2741

Product Information. U.S. and Canada: 1.800.992.1776

Outside U.S. and Canada: www.metro.com/contactus

LO2-243

Printed in U.S.A. 3/20

Information and specifications are subject to change without notice. Please confirm at time of order.

Copyright © 2020 InterMetro Industries Corp.