

# Certificate of Quality



|                          |                                       |
|--------------------------|---------------------------------------|
| <b>Part Number</b>       | <b>MFLX06411-03</b>                   |
| <b>Lot Number</b>        | <b>28695756</b>                       |
| <b>Part Description</b>  | <b>TUBING SILICONE 3/16X5/16" 25'</b> |
| <b>Country of Origin</b> | <b>USA</b>                            |

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: Peroxide Cured Silicone
- Per the manufacturers technical data sheet:
- FDA Approved for food contact
- USP Class VI - OK
- Autoclavable - OK
- This product is not approved or intended for, and should not be used for medical, clinical, surgical or other patient oriented applications.
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.