



Material	48372-540
Material description	Iodine (Iodine-Iodide) 0.1N
Grade	Volumetric Solution
Molecular formula	I ₂
Molecular mass	253,81
Additional information	

Characteristics

Specifications

Normality	0,0995 --> 0,1005 N
Appearance	Passes test
Traceable to NIST SRM 136f	Confirmed

Signature

Additional information

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada

Transparent, reddish brown solution, no residue, precipitate, floater