



**VWR® GENERAL-PURPOSE LABORATORY LABELING TAPE**

As per the manufacturer, the below product meets the following criteria:

<b>North American Catalog No:</b>	89097-930
<b>Lot Number:</b>	706796
<b>Description:</b>	VWR Tape Grn 1/2x500 Pk6rl
<b>Date of Manufacture:</b>	June 22, 2020
<b>Product Shelf life:</b>	18 months
<b>Country of Origin:</b>	USA

**QC Testing**

The following data are nominal values based on PSTC, ASTM and other standard tests. The data should not be considered specifications.

Adhesion to Steel (oz/in of width) 50 (13.7 N/25mm) PSTC-101, Quick Stick (oz/in of width) To Kraft 8 (2.2 N/25mm), To Steel 23 (6.3 N/25mm)- Tensile Strength (lbs/in of width) 31 (136 N/25mm) PSTC- 131- Elongation (% at break) 3 PSTC-131

-Total Thickness (mils) 6.3 (0.16mm)

**Product Specifications**

**Material:** Flatback paper with synthetic rubber adhesive

**BSE/TSE:** The above product does not contain materials infected with or containing either BSE (bovine spongiform encephalopathy) or TSE (transmittable spongiform encephalopathy) a.k.a. Mad Cow Disease.

**Non-Pyrogenic:** N/A

**DNase & RNase Free:** NA

**Cytotoxicity:** N/A

**Latex Free:** The above listed does not contain raw natural rubber latex (NRL) or dried natural rubber (DNR).

**Signed:**

Jamie Ethier  
VP Global Quality  
VWR, Part of Avantor

**Date:** October 14, 2020