


SMZ171 Series

SMZ-171-BP





SMZ161 Series

SMZ-171-BP

Optical system	Greenough
Observation tube	Binocular head
Inclination	45° inclined
Interpupillary distance	52-75mm
Diopter adjustment	On both eyepieces, +/- 5 diopter
Eyepieces	Widefield N-WF10X/23mm with diopter adjustment
Objectives system	Zoom, ratio 6.7:1
Objectives	0.75X-5X
Working distance	110mm
Stand type	Pole type stand without illumination
Head holder	For Ø32mm column and Ø76mm head
Focus mechanism	Coarse focusing system with tension adjustment
Focusing stroke	49mm
Accessories included	Dust cover, black/white stage plate
Dimensions LxWxH	303x239x370mm
Base	303x239mm
Column/Arm	Height 290mm, Ø32mm
Net weight	5.45kg