



Sodium Hydroxide, 1N

Volumetric Solution

BAKER ANALYZED Reagent

Product No. 5635
Specifications current as of:
Jun 8 2011

TEST	SPECIFICATION
Standardization at 25°C traceable to NIST Standard Reference Material.	
SRM No...	
Appearance (Clear, colorless liquid)	Passes Test
Normality (eq/L)	0.995 - 1.005
Trace Impurities (in ppm):	
Chloride (Cl)	max. 5
Heavy Metals (as Pb)	max. 1
Iron (Fe)	max. 0.5
Filtered through a 0.2 micron filter.	
For Laboratory, Research or Manufacturing Use	
Protect from air to avoid absorption of carbon dioxide.	