



VWR® Centrifuge Tubes with Flat or Plug Caps

Description: Centrifuge Tubes with Flat or Plug Caps, Sterilized/Non-Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.


Materials

Tube: PP (Polypropylene), Compliance with USP Class VI Color: Clear

Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: White

Rack: Paper Color: White

Dimensions (Unit:mm)



Features

- Easy-to-read black graduations in $\pm 2\%$ increments
- Contains a large, white frosted writing area
- Autoclavable at 121°C and freezable to -80°C
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic
- Meet the requirements of RoHS
- Leak-proof
- DNase/RNase free
- Sterilized by irradiation SAL10⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® Centrifuge Tubes with Flat or Plug Caps

VWR EU EAN, No.	VWR NA Cat. No.	Capacity (ml)	Maximum RCF	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
525-0604	10025-686	15ml	12,500g	Conical-Bottom	Clear	Flat cap	White	Yes	50/Bag, 10 Bags/CS, 500/CS
525-0605	10025-688	15ml	12,500g	Conical-Bottom	Clear	Plug cap	White	Yes	50/Bag, 10 Bags/CS, 500/CS
525-0606	10025-690	15ml	12,500g	Conical-Bottom	Clear	Flat cap	White	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-0607	10025-692	15ml	12,500g	Conical-Bottom	Clear	Plug cap	White	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-0608	10025-694	50ml	12,500g	Conical-Bottom	Clear	Flat cap	White	Yes	25/RK/PK, 20 PK/CS, 500/CS
525-0609	10025-696	50ml	12,500g	Conical-Bottom	Clear	Plug cap	White	Yes	25/RK/PK, 20 PK/CS, 500/CS
525-0610	10025-698	50ml	12,500g	Conical-Bottom	Clear	Flat cap	White	Yes	50/Bag, 10 Bags/CS, 500/CS
525-0612	10025-702	50ml	12,500g	Conical-Bottom	Clear	Plug cap	White	Yes	50/Bag, 10 Bags/CS, 500/CS
525-0611	10025-700	50ml	12,500g	Freestanding	Clear	Flat cap	White	Yes	50/Bag, 10 Bags/CS, 500/CS
525-0633	10025-712	50ml	10,000g	Freestanding	Clear	Flat cap	White	No	50/Bag, 10 Bags/CS, 500/CS
525-0634	10026-072	50ml	10,000g	Freestanding	Clear	Flat cap	White	Yes	50/Bag, 10 Bags/CS, 500/CS
525-0635	10026-074	50ml	10,000g	Freestanding	Clear	Flat cap	White	No	Bulk, 500/CS
	75845-202	50ml	12,500g	Conical-Bottom	Clear	Flat cap	White	Yes	Tubes:50/bag, Caps:500/Bag, 500/CS
	75845-204	50ml	12,500g	Freestanding	Clear	Flat cap	White	Yes	Tubes:50/bag, Caps:500/Bag, 500/CS
	75845-206	50ml	12,500g	Conical-Bottom	Clear	Plug cap	White	Yes	Tubes:50/bag, Caps:500/Bag, 500/CS
525-0636	10026-076	15ml	7,500g	Conical-Bottom	Clear	Flat cap	White	Yes	50/Bag, 10 Bags/CS, 500/CS
525-0637	10026-078	50ml	7,500g	Conical-Bottom	Clear	Flat cap	White	Yes	50/Bag, 10 Bags/CS, 500/CS
	76026-480	15ml	12,500g	Conical-Bottom	Clear	Flat cap	White	No	Bulk, 500/CS
	76026-482	15ml				Flat cap	Blue	No	500/Bag, 10 Bags/CS, 5000/CS
	76319-832	15ml				Flat cap	White	Yes	500/Bag, 10 Bags/CS, 5000/CS
	76319-834	15ml				Plug cap	White	Yes	500/Bag, 10 Bags/CS, 5000/CS
	76319-836	50ml				Flat cap	White	Yes	500/Bag, 10 Bags/CS, 5000/CS
	76319-830	50ml				Plug cap	White	Yes	500/Bag, 10 Bags/CS, 5000/CS