

## Certificate Of Analysis

**Whatman™**

### Product Information

Grade: 1CHR

Product Number: 3001-861

Product Name: GR 1CHR 20x20CM 100/PK

Lot Number/Serial Number: 17208758

Date of Manufacture: 2-Nov-2020

### Quality Control

| Test                    | Range   | Results |
|-------------------------|---------|---------|
| Grammage,gsm            | 80-93   | 84      |
| Thickness, um/53kpa     | 160-196 | 167     |
| Porosity, S/200ml/sq.in | 13-31   | 19      |

### Conformance & Quality systems statement

This is to certify that this product conforms to GE Healthcare Life Sciences specifications.

All products are manufactured via a Quality System certified to ISO9001:2015 and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

### Electronic signature

This document has electronically produced and is valid without a signature.

Version AD



imagination at work

[www.gelifesciences.com](http://www.gelifesciences.com)

Global Life Sciences Solutions Operations UK Ltd.

Amersham Place Little Chalfont

Buckinghamshire HP7 9NA UK