

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® 11 mm Snap Ring Vials, Kits, Caps, and Inserts, Wide Opening
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2023-08-11

Vials

VWR Catalogue Number	Description	Capacity	Bottom style	Material	Color
82028-412A	Vials with insert	0.30 ml	Conical	PP	Natural
10803-890A	Vials with fused insert	0.35 ml	Conical	Borosilicate glass	Clear
66009-822A	Vials	2.00 ml	Flat bottom	Borosilicate glass	Clear
82028-414A	Vials with ID patch	2.00 ml	Flat bottom	Borosilicate glass	Clear
82028-416A	Vials with ID patch	2.00 ml	Flat bottom	Borosilicate glass	Amber

Clear Vial kits with 2.00 ml capacity and Flat bottom

VWR Catalogue Number	Material	Septum	Cap included	Cap size	Color material	Cap color
66009-868A	Borosilicate glass	Clear PTFE/synthetic red rubber	yes	11 mm	PE	Natural
66009-872A	Borosilicate glass	Red PTFE/white silicone	yes	11 mm	PE	Natural

Caps - Snap caps with septa

VWR Catalogue Number	Material	Septum	Cap color	Cap size
46610-718A	PE	Clear PTFE/Synthetic red rubber	Natural	11 mm
89239-010A	PE	Clear PTFE/Synthetic red rubber	Blue	11 mm
46610-720A	PE	Red PTFE/White silicone	Natural	11 mm

89239-012A	PE	Red PTFE/White silicone	Blue	11 mm
82028-418A	PE	Blue PTFE/White silicone, pre-slit	Natural	11 mm
89239-014A	PE	Blue PTFE/White silicone, pre-slit	Blue	11 mm

Clear Inserts without cap

VWR Catalogue Number	Description	Capacity	Bottom style	Material
82028-450A	Insert with poly-support spring for wide opening vials	0.33 ml	Conical bottom	PP
46610-762A	Insert with poly-support spring for wide opening vials	0.38 ml	Conical bottom	Glass
82028-452A	Insert, pulled point bottom, for wide opening vials	0.40 ml	Conical bottom	Glass
82028-454A	Insert for wide opening vials	0.50 ml	Flat bottom	Glass

Quality System Compliance

Products are manufactured under the [ISO 9001:2015](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Refer Table.
Sterilization	Not Applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not Applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not Applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Not Applicable.
Shelf life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	August 2023