



| | |
|----------------------|-----------------|
| Material | 67482-690 |
| Material description | Petroleum ether |
| Grade | ACS |
| CAS Number | 8032-32-4 |

Characteristics

Specifications

| | |
|---|--------------|
| Acidity | Passes test |
| Boiling range (35-60°C) | Passes test |
| Color (Pt-Co) | Max. 10 APHA |
| Residue after evaporation | Max. 0.001 % |
| Meets MIL-PRF-680 Type 1 specifications | Conforms |

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.