

F32036 0,1M Potassium bromideFIX1,00SA

00000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

SECTION 1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Potassium bromide concentrate

Number : 00000021505

Product Use Description : Laboratory chemicals

Manufacturer or supplier's details : Honeywell International Inc.
1953 South Harvey Street
Muskegon, MI 49442

For more information call : 1-800-368-0050
+1-231-726-3171

(Monday-Friday, 9:00am-5:00pm)

In case of emergency call : **Medical: 1-800-498-5701 or +1-303-389-1414**
: **Transportation (CHEMTREC): 1-800-424-9300 or**
: **+1-703-527-3887**
:
: (24 hours/day, 7 days/week)

SECTION 2. HAZARDS IDENTIFICATION**Emergency Overview**

Form : liquid

Color : colourless

Odor : odourless

Classification of the substance or mixture**|| Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).****|| Precautionary statements : Prevention:**
Wear protective gloves/protective clothing/eye protection/face protection.

F32036 0,1M Potassium bromideFIX1,00SA

000000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

Hazards not otherwise classified : Repeated or prolonged exposure may irritate eyes, skin and respiratory system.

Carcinogenicity

No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP, IARC, or OSHA.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Formula : BrK

Chemical nature : Mixture

Chemical name	CAS-No.	Concentration
Water	7732-18-5	>=80.00 - <=90.00 %
Potassium bromide	7758-02-3	>=10.00 - <=20.00 %

SECTION 4. FIRST AID MEASURES

General advice : First aider needs to protect himself. Move out of dangerous area. Take off all contaminated clothing immediately.

Inhalation : Remove to fresh air. If respiratory problems develop, obtain medical attention.

Skin contact : Wash off with soap and plenty of water. If skin irritation persists, call a physician.

Eye contact : Rinse immediately with plenty of water, also under the eyelids. If eye irritation persists, consult a specialist.

Ingestion : Rinse mouth with water. Do not induce vomiting without medical advice. Never give anything by mouth to an unconscious person.

F32036 0,1M Potassium bromideFIX1,00SA

00000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

SECTION 5. FIREFIGHTING MEASURES

- Suitable extinguishing media : Water spray
Foam
Carbon dioxide (CO₂)
Dry powder
- Specific hazards during firefighting : In case of fire hazardous decomposition products may be produced such as:
Hydrogen bromide
Bromine compounds
Potassium oxide
- Special protective equipment for firefighters : In the event of fire, wear self-contained breathing apparatus.
- Further information : Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

SECTION 6. ACCIDENTAL RELEASE MEASURES

- Personal precautions : Evacuate personnel to safe areas.
Wear personal protective equipment. Unprotected persons must be kept away.
Ensure adequate ventilation.
Avoid contact with skin, eyes and clothing.
- Environmental precautions : Prevent further leakage or spillage if safe to do so.
Do not let product enter drains.
If the product contaminates rivers and lakes or drains inform respective authorities.
- Methods for cleaning up : Soak up with inert absorbent material.
Sweep up and shovel into suitable containers for disposal.
Dispose of in accordance with local regulations.

SECTION 7. HANDLING AND STORAGE**Handling**

- Handling : Wear personal protective equipment.
Use only in well-ventilated areas.
- Advice on protection against : Normal measures for preventive fire protection.

F32036 0,1M Potassium bromideFIX1,00SA

00000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

fire and explosion

Storage

Requirements for storage areas and containers : Store in original container.
Keep containers tightly closed in a dry, cool and well-ventilated place.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Protective measures : Ensure that eyewash stations and safety showers are close to the workstation location.
Legal requirements are to be considered in regard of the selection, use and care of personal protective equipment.
Do not breathe vapours or spray mist.
Avoid contact with skin, eyes and clothing.

Engineering measures : Ensure adequate ventilation, especially in confined areas.

Eye protection : Safety glasses with side-shields

Hand protection : Latex gloves
Gloves must be inspected prior to use.
Replace when worn.

Skin and body protection : Wear suitable protective equipment.

Respiratory protection : In case of insufficient ventilation wear suitable respiratory equipment.

Hygiene measures : General industrial hygiene practice.

Exposure Guidelines

Contains no substances with occupational exposure limit values.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Physical state : liquid

Color : colourless

Odor : odourless

pH : Note: no data available

SAFETY DATA SHEET



F32036 0,1M Potassium bromideFIX1,00SA

00000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

Melting point/range	: Note: no data available
Boiling point/boiling range	: Note: no data available
Flash point	: Note: Not applicable
Lower explosion limit	: Note: Not applicable
Upper explosion limit	: Note: Not applicable
Density	: 1.0 g/cm ³ at 20 °C
Water solubility	: Note: completely miscible
Partition coefficient: n-octanol/water	: Note: no data available
Decomposition temperature	: Note: No decomposition if stored and applied as directed.
Viscosity, dynamic	: Note: no data available
Oxidizing properties	: The substance or mixture is not classified as oxidizing.
Molecular weight	: 119.00 g/mol

SECTION 10. STABILITY AND REACTIVITY

Chemical stability	: Stable under recommended storage conditions.
Possibility of hazardous reactions	: No dangerous reaction known under conditions of normal use.
Conditions to avoid	: None known.

F32036 0,1M Potassium bromideFIX1,00SA

00000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

Incompatible materials to avoid : Strong oxidizing agents
Strong acids

Hazardous decomposition products : In case of fire hazardous decomposition products may be produced such as:
Hydrogen bromide
Bromine compounds
Potassium oxide

SECTION 11. TOXICOLOGICAL INFORMATION

Acute oral toxicity : Note: no data available

Acute inhalation toxicity : Note: no data available

Acute dermal toxicity : Note: no data available

Skin irritation : Note: no data available

Eye irritation : Note: no data available

Sensitisation : Note: no data available

SECTION 12. ECOLOGICAL INFORMATION**Ecotoxicity effects**

Toxicity to fish : Note: no data available

Toxicity to daphnia and other aquatic invertebrates : Note: no data available

Toxicity to algae : Note: no data available

F32036 0,1M Potassium bromideFIX1,00SA

00000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

Elimination information (persistence and degradability)

Biodegradability : Note: The methods for determining biodegradability are not applicable to inorganic substances.

Further information on ecology**SECTION 13. DISPOSAL CONSIDERATIONS**

Disposal methods : Observe all Federal, State, and Local Environmental regulations.

SECTION 14. TRANSPORT INFORMATION

DOT Not dangerous goods

TDG Not dangerous goods

IATA Not dangerous goods

IMDG Not dangerous goods

SECTION 15. REGULATORY INFORMATION**Inventories**

US. Toxic Substances Control Act : On TSCA Inventory

Australia. Industrial Chemical (Notification and Assessment) Act : On the inventory, or in compliance with the inventory

Canada. Canadian Environmental Protection Act (CEPA). Domestic Substances List (DSL) : All components of this product are on the Canadian DSL

Japan. Kashin-Hou Law List : On the inventory, or in compliance with the inventory

Korea. Existing Chemicals Inventory (KECI) : On the inventory, or in compliance with the inventory

Philippines. The Toxic : On the inventory, or in compliance with the inventory

F32036 0,1M Potassium bromideFIX1,00SA

000000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

Substances and Hazardous
and Nuclear Waste Control
Act

China. Inventory of Existing Chemical Substances : On the inventory, or in compliance with the inventory

New Zealand. Inventory of Chemicals (NZIoC), as published by ERMA New Zealand : On the inventory, or in compliance with the inventory

National regulatory information

SARA 302 Components : No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components : This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards : Acute Health Hazard

California Prop. 65 : This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

New Jersey RTK : Potassium bromide 7758-02-3

Pennsylvania RTK : Potassium bromide 7758-02-3

SECTION 16. OTHER INFORMATION

	HMIS III	NFPA
Health hazard	: 1	1
Flammability	: 0	0
Physical Hazard	: 0	
Instability	:	0

F32036 0,1M Potassium bromideFIX1,00SA

00000021505

Version 2.0

Revision Date 07/31/2017

Print Date 10/17/2019

Hazard rating and rating systems (e.g. HMIS® III, NFPA): This information is intended solely for the use of individuals trained in the particular system.

Further information

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text. Final determination of suitability of any material is the sole responsibility of the user. This information should not constitute a guarantee for any specific product properties.

Changes since the last version are highlighted in the margin. This version replaces all previous versions.

Prepared by Honeywell Performance Materials and Technologies Product Stewardship Group