



150PR

G-Biosciences ♦ 1-800-628-7730 ♦ 1-314-991-6034 ♦ technical@GBiosciences.com

A Geno Technology, Inc. (USA) brand name


OneQuant™ HOOK™ Sulfo- NHS-Biotin

(Cat.# 786-698)



think proteins! think G-Biosciences www.GBiosciences.com

INTRODUCTION



Sulfo-NHS-Biotin is a commonly used amine reactive biotin-labeling reagent. OneQuant™ HOOK™ Sulfo-NHS-Biotin are single aliquots of Sulfo-NHS-Biotin reagent in convenient, premeasured microtubes. Each microtube contains 1mg of the biotin reagent. Simply dissolve one vial with appropriate solvent and add to the reaction.

Formula: C₁₄H₁₈O₈N₃S₂Na
Molecular Weight: 443.43
Spacer Arm Length (Å): 13.5
Reaction Group: sulfo-NHS-ester
Membrane Permeable: NO
Water Soluble: YES
Cleavable/ Reversible: NO
Reaction pH: 7-9
CAS No: 119616-38-5

ITEMS SUPPLIED

Description	Size
OneQuant™ HOOK™ Sulfo-NHS-Biotin	1mg/vial, 1 strip of 8 vials

STORAGE CONDITIONS

It is shipped at ambient temperature. Upon receipt store product at 4°C protected from moisture. Stable for one year.

RELATED PRODUCTS

Download our Protein Labeling & Conjugation and Antibody Production Handbook.



<http://info2.gbiosciences.com/complete-protein-labeling-conjugation-handbook>

<http://info2.gbiosciences.com/complete-antibody-production-handbook>

For other related products, visit our website at www.GBiosciences.com or contact us.

Last saved: 7/14/2015 TNN



www.GBiosciences.com