

Material No.: 7581-60

:

Revision No.: 0

Product Specification

Meets ACS Reagent Chemical Requirements,

Test	Specification
ACS - Assay (NaCl) (by Ag titrn)	≥ 99.0 %
ACS - pH of 5% Solution at 25°C	5.0 - 9.0
ACS - Insoluble Matter	≤ 0.005 %
ACS - Iodide (I)	≤ 0.002 %
ACS - Bromide (Br)	≤ 0.01 %
ACS - Chlorate and Nitrate (as NO ₃)	≤ 0.003 %
ACS - Phosphate (PO ₄)	≤ 5 ppm
ACS - Sulfate (SO ₄)	≤ 0.004 %
ACS - Barium (Ba)	Passes Test
ACS - Heavy Metals (as Pb)	≤ 5 ppm
ACS - Iron (Fe)	≤ 2 ppm
ACS - Calcium (Ca)	≤ 0.002 %
ACS - Magnesium (Mg)	≤ 0.001 %
ACS - Potassium (K)	≤ 0.005 %

For Laboratory, Research, or Manufacturing Use
Appearance (white, crystalline granules)

Packaging Site: Paris Mfg Ctr & DC