



Acetic Acid, Glacial

USP, FCC

Product No. 9522
 Specifications current as of: Jun 8 2011

TEST	SPECIFICATION
GMP Manufactured Product	
Meets U.S.P. & F.C.C. Requirements	
Identification	Passes Test
Congealing Temperature (solidification point)	min. 15.6 °C.
Nonvolatile Residue (USP), mg/20 mL	max. 1.0
Nonvolatile Residue (FCC)	max. 0.005 %
Chloride (Cl)	Passes Test
Sulfate (SO ₄)	Passes Test
Heavy Metals (as Pb)	max. 5 ppm
Readily Oxidizable Substances	Passes Test
Assay (CH ₃ COOH)	99.5 - 100.5 %
Lead (Pb)	max 0.5 mg/kg
Only Class 3 Solvents are likely to be present. Residual concentration of all Class 3 Solvents is <0.35%.	
Bulk Food Chemical	
Bulk Pharmaceutical Chemical	
CAUTION: For manufacturing, processing or repacking.	
For a Certificate of Lot Analysis, visit www.avantormaterials.com .	
IMPORTANT: Material will freeze if stored below 17°C (63°F).	
Preserve in tight containers, and store at room temperature	