

Thermo Scientific GC consumables for Agilent GC & GC-MS Instruments

GC & GC/MS: 5890, 6890, 7820, 7890

Injection Port Liners



Description	Agilent Part Number	VWR equivalent Catalog No. (5 Pack)	VWR equivalent Catalog No. (25 Pack)
Split/Splitless FocusLiner with Double Taper, 4 x 6.3 x 78.5mm	5183-4709	10055-344	10171-678
Split Straight Liner 4 x 6.3 x 78.5mm, Quartz Wool	5183-4691	10055-362	10055-332
Split/Splitless FocusLiner 4 x 6.3 x 78.5mm, Quartz Wool	210-4004-5	10055-330	10055-328
Split/Splitless FocusLiner with Single Taper, 4 x 6.3 x 78.5mm	5183-4693	10055-342	10055-340
Split/Splitless Liner with Single Taper 4 x 6.3 x 78.5mm	5183-4695	10055-346	10055-376
Split/Splitless FocusLiner with Double Taper, 4 x 6.3 x 78.5mm	5183-3315	10055-348	10055-378

Ferrules



Description	Agilent Part Number	VWR equivalent Catalog No.	Pack Size
100% Graphite Ferrule for 0.1-0.32mm ID Columns	500-2114 / 5080-8853	10055-910	10
100% Graphite Ferrule for 0.45-0.53mm ID Columns	500-2118	10055-912	10
15% Graphite/85% Vespel Ferrules for 0.1-0.25mm ID Columns	5181-3323	10055-746	10
15% Graphite/85% Vespel Ferrules for 0.32mm ID Columns	5062-3514	10055-748	10
SilTite Kit for 0.1-0.25mm ID Columns (GC/MS Interface)	5188-5361	10055-082	1
SilTite Kit for 0.32mm ID Columns (GC/MS Interface)	5188-5362	10055-084	1

Septa



Description	Agilent Part Number	VWR equivalent Catalog No.	Pack Size
BTO Low Bleed Silicone Septa 11mm Diameter	5183-4757	10141-486	50
TR-Green Low Coring Silicone Septa 11mm Diameter	5183-4759	10056-088	50

Syringes



Description	Agilent Part Number	VWR equivalent Catalog No.
5µL Fixed Needle Syringe, 42mm Length, 23 Gauge, Cone Needle	9301-0892	10055-192
5µL Fixed Needle Syringe, 42mm Length, 23-26 Gauge, Cone Needle	5182-0835	10055-196
10µL Fixed Needle Syringe, 42mm Length, 23 Gauge, Cone Needle	9301-0713	10055-212
10µL Fixed Needle Syringe, 42mm Length, 23-26 Gauge, Cone Needle	5181-1267	10055-218

Inlet Base Seals



Description	Agilent Part Number	VWR equivalent Catalog No. (10 Pack)	VWR equivalent Catalog No. (2 Pack)
Gold Base Seals 0.8mm ID	5190-2209	10055-762	10055-764
Siltek Treated Base Seals 0.8mm ID	–	10055-770	10055-772

TraceGOLD GC Columns



Description	Agilent Part Number	VWR equivalent Catalog No.
TG-5MS 30m × 0.25mm × 0.25µm	122-5532	10038-950
TG-5MS 60m × 0.25mm × 0.25µm	122-5562	10038-956
TG-5MS 30m × 0.32mm × 0.25µm	123-5532	10038-966
TG-1MS 30m × 0.25mm × 0.25µm	122-0132	10056-756
TG-1MS 60m × 0.25mm × 0.25µm	122-0162	10056-762
TG-WaxMS 30m × 0.25mm × 0.25µm	122-7032	10055-728
TG-WaxMS 60m × 0.25mm × 0.25µm	122-7062	10056-392
TG-5HT 30m × 0.25mm × 0.25µm	123-7032	10056-600

SWAP-IT GC-MS Interface



Description	VWR Catalog No.
SWAP-IT GC-MS Interface Kit	10055-994
SWAP-IT GC-MS Interface – SWAP-IT Only	10055-996

Gas Filters



Description	Agilent Part Number	VWR equivalent Catalog No.
Super Clean Moisture Cartridge Filter	CP17971	10056-206
Super Clean Oxygen Cartridge Filter	CP17970	10056-266

Equipment



Description	Agilent Part Number	VWR equivalent Catalog No.
GFM Pro Gas Flow Meter	5067-0223	10057-406
GLD Pro Gas Leak Detector	G3388A	10057-298
Tool Kit for Agilent Instruments	5180-4162	10055-138

Vials & Closures



Description	VWR equivalent Catalog No.
AVCS Kit, Clear Screw Thread Vials, Red Rubber/Clear PTFE Seals 100pk	89239-444
AVCS Kit, Amber Screw Thread Vials, Red Rubber/Clear PTFE Seals 100pk	10119-726

VWR  **1.800.932.5000**
We Enable Science **vwr.com**



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2017 VWR International, LLC. All rights reserved.

0117 Lit. No. 201942W