

VWR, part of Avantor
100 Matsonford Road
Building 1, Suite 200
Radnor PA 19087
t: 1.800.932.5000
www.vwr.com

CUBITAINER 20L (5 Gal)

- Description:** A Cubitainer is a plastic liner inside a corrugated fiberboard box. The liner acts as a bladder, collapsing as liquid is dispensed.
- Material:** LDPE liner with corrugated fiberboard overpack and LDPE tap included.
- Certification:** USA/+AA9711
4G/Y32.6/S/18
- Closure:** 32 mm threaded opening with LDPE dispensing tap.
Product is shipped with a cap and a detached spigot (or tap). The customer must remove the cap to attach the tap.
- Dimensions:** 12 ¼ x 12 ¼ x 12 ¼ inches
- Images:** Inner liner and assembled inner liner with outer overpack.



Tap (included)



For more information, please refer to the specific product page on vwr.com or contact:

VWR Technical Product Support
TechnicalProductSupportNA@vwr.com
Tel: 888.897.5463