

# AmpliPur™

EXPERT TIPS



EXCELLENT FIT FILTER TIPS FOR MOLECULAR BIOLOGY APPLICATIONS

| Excellent Fit

| Ideal for Molecular Biology Applications

| ISO 8655 Compliant




 **VWR™**  
part of avantor

order on  
[VWR.COM](http://VWR.COM)



## Excellent Fit\*

AmpliPur Expert Tips are designed to fit various pipette brands for maximum convenience on the bench. Equipped with a filter, they also help prevent cross-contamination.



## Ideal for Molecular Biology Applications

These tips significantly contribute to molecular testing performance by ensuring sample integrity during the pipetting process.

## ISO 8655 Compliant

Manufactured according to the certified Gilson quality standard, AmpliPur Expert Tips comply with ISO 8655 requirements and are certified free of RNase and DNase, preserving your sample integrity.



\*Pipettes compatibility need to be verified by the user.

| Technical Information |  |
|-----------------------|--|
| Gilson Quality System | ISO 9001   |
| Product Compliance    | ISO 8655   |
| Analysis Certificate  | Free of RNase and DNase<br>Hazardous Materials (2015/863/EU) |
| Materials             | <b>Tips:</b> Polypropylene (PP)                              |
|                       | <b>Filter:</b> Polyethylene (PE)                             |
|                       | <b>Rack:</b> Polypropylene (PP)                              |
|                       | <b>Heat shrink:</b> Plastic Optical Fibre (POF)              |
| Packaging             | 10 racks of 96 tips (960 tips)                               |
| Disposable            | <b>Tips:</b> Yes   |
|                       | <b>Rack:</b> Reusable  |
| Warranty              | 3 years  |

| AmpliPur™ Expert Tips |   |                                |
|-----------------------|---|--------------------------------|
| Volume Range          | Nominal Volume of Compatible Pipettes** | Part Number (Commercial Units) |
| 0.5-10 µL             | 2 µL / 10 µL                            | 76514-164                      |
| 1-20 µL               | 20 µL                                   | 76514-166                      |
| 100-1000 µL           | 1000 µL                                 | 76514-168                      |
| 100-1200 µL           | 1000 µL / 1200 µL                       | 76514-170                      |

\*\*Within ISO 8655 specifications.

\*\*\*Compatible with many of 100 µL pipettes, except Gilson P100 single channel.

AmpliPur Expert Tips are only available in select countries. Please contact your local Gilson representative for more information.

