

# Certificate of Quality

<b>VWR Catalogue Number</b>	See List
<b>Description</b>	VWR® Layflat Poly Bags
<b>Country of Origin</b>	Manufactured in USA
<b>Date of Issue (yyyy-mm-dd)</b>	2021-04-20

VWR Catalogue Number	WxH	Thickness [mil]
89071-124	12.7 x 20.3 cm (5 x 8")	1
89071-228	30.5 x 30.5 cm (12 x 12")	1
89071-240	30.5 x 38.1 cm (12 x 15")	1
89071-264	35.6 x 50.8 cm (14 x 20")	1
89071-290	45.8 x 61 cm (18 x 24")	1
89071-306	61 x 61 cm (24 x 24")	1
89071-296	50.8 x 61 cm (20 x 24")	1
89071-492	20.3 x 30.5 cm (8 x 12")	1.25
89069-122	12.7 x 17.7 cm (5 x 7")	1.5
89069-150	15.2 x 35.6 cm (6 x 14")	1.5
89069-160	17.7 x 20.3 cm (7 x 8")	1.5
89069-204	20.3 x 35.6 cm (8 x 14")	1.5
89069-210	20.3 x 61 cm (8 x 24")	1.5
89069-232	25.4 x 30.5 cm (10 x 12")	1.5
89069-252	30.5 x 30.5 cm (12 x 12")	1.5
89069-262	30.5 x 35.6 cm (12 x 14")	1.5
89069-068	5.2 x 17.7 cm (2 x 7")	1.5
89069-488	10.2 x 10.2 cm (4 x 4")	2
89069-500	10.2 x 35.6 cm (4 x 14")	2
89069-542	15.2 x 20.3 cm (6 x 8")	2
89069-556	15.2 x 30.5 cm (6 x 12")	2
89069-596	20.3 x 30.5 cm (8 x 12")	2

89069-614	22.9 x 30.5 cm (9 x 12")	2
89069-668	30.5 x 40.6 cm (12 x 16")	2
89069-686	30.5 x 61 cm (12 x 24")	2
89069-752	50.8 x 76.2 cm (20 x 30")	2
89069-768	66 x 91.4 cm (26 x 36")	2
89069-788	71.1 x 91.4 cm (28 x 36")	2
89069-814	81.2 x 122 cm (32 x 48")	2
89069-826	91.4 x 152.4 cm (36 x 60")	2
89069-930	15.2 x 20.3 cm (6 x 8")	3
89069-986	22.9 x 40.6 cm (9 x 16")	3
89070-034	30.5 x 45.8 cm (12 x 18")	3
89070-038	30.5 x 61 cm (12 x 24")	3
89070-070	40.6 x 61 cm (16 x 24")	3
89070-084	45.8 x 61 cm (18 x 24")	3
89070-104	61 x 61 cm (24 x 24")	3
89070-126	71.1 x 91.4 cm (28 x 36")	3
89070-130	76.2 x 91.4 cm (30 x 36")	3
89070-518	101.6 x 106.7 cm (40 x 42")	4
89070-242	12.7 x 12.7 cm (5 x 5")	4
89070-252	12.7 x 20.3 cm (5 x 8")	4
89070-286	15.2 x 30.5 cm (6 x 12")	4
89070-292	15.2 x 61 cm (6 x 24")	4
89070-308	20.3 x 25.4 cm (8 x 10")	4
89070-314	20.3 x 30.5 cm (8 x 12")	4
89070-316	20.3 x 35.6 cm (8 x 14")	4
89070-312	20.3 x 71.1 cm (8 x 28")	4
89070-350	27.9 x 40.6 cm (11 x 16")	4
89070-368	30.5 x 38.1 cm (12 x 15")	4
89070-370	30.5 x 40.6 cm (12 x 16")	4
89070-376	30.5 x 61 cm (12 x 24")	4
89070-406	40.6 x 45.8 cm (16 x 18")	4
89070-180	5.2 x 7.6 cm (2 x 3")	4

89070-444	50.8 x 61 cm (20 x 24")	4
89070-446	50.8 x 76.2 cm (20 x 30")	4
89070-462	61 x 76.2 cm (24 x 30")	4
89070-470	61 x 91.4 cm (24 x 36")	4
89070-200	7.6 x 10.2 cm (3 x 4")	4
89070-202	7.6 x 12.7 cm (3 x 5")	4
89070-492	76.2 x 122 cm (30 x 48")	4
89070-488	76.2 x 91.4 cm (30 x 36")	4
89070-514	91.4 x 122 cm (36 x 48")	4
89072-482	61 x 152.4 cm (24 x 60")	6
89072-476	61 x 91.4 cm (24 x 36")	6

### Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

### QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

### Product Specifications

<b>Material</b>	Low Density Polyethylene (LDPE), Linear Low-Density Polyethylene (LLDPE).
<b>Sterilization</b>	Not applicable.
<b>ATP Assay</b>	Not applicable.
<b>DNase &amp; RNase Free</b>	This product is free of any detectable DNase/RNase contamination.
<b>BSE/TSE</b>	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
<b>Cytotoxicity</b>	Not applicable.
<b>Non-Pyrogenic Statement</b>	Not applicable.

<b>Latex Statement</b>	The product is Latex free.
<b>BPA Statement</b>	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
<b>DEHP Statement</b>	Please see RoHS information below.
<b>RoHS</b>	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr <sup>6+</sup> ), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP), DEHP (di(2- ethylhexyl) phthalate)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
<b>REACH Statement</b>	No substances as listed on >Substance of Very High Concern< candidate list published 01/19/2021 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
<b>Storage Conditions</b>	Recommended indoor storage condition / out of elements. Temperature control is not necessary.
<b>Shelf Life</b>	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.