

Product Specification

Catalog Number: BDH9268

Catalog Description: BDH POTASS PHOSPHAT MONO

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics

Specifications

Chloride	<= 0.001 %
Heavy Metals	<= 0.001 %
Insolubles	<= 0.01 %
Iron	<= 0.002 %
Loss on Drying	<= 0.2 %
Nitrogen Compounds	<= 0.001 %
pH (5%, Water) @25C	4.1 - 4.5
Purity	>= 99.0 %
Sodium	<= 0.005 %
Sulfate	<= 0.003 %

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager