

CERTIFICATE OF STERILITY



Veltek Associates, Inc.

LABORATORIES

PRODUCT STERILITY DETAILS

The following lot numbers have passed the current USP Test for Sterility.

Date: **25-SEP-2019**

Reference #: **V12950**

Product Code: **DSTER-WFI-B-70**

Description: **STER-AHOL 70% WFI Gallon**

PRODUCT LOT

Lot #19-DSTER-910874	EXP	01-AUG-2022
	EXP	
	EXP	
	EXP	
	EXP	
	EXP	
	EXP	

IRRADIATION#2505535
106 CASES

A handwritten signature in black ink, appearing to read "M. Lang".

25 Sep 2019

Authorized Signature-Quality Assurance

Date

Veltek Associates, Inc.

15 Lee Blvd • Malvern, PA 19355-1234
Tel: (610) 644-8335 Fax: (610) 644-8336

CERTIFICATE OF IRRADIATION



Veltek Associates, Inc.

PRODUCT IRRADIATION DETAILS

Irradiation Date: **01-SEP-2019**

Reference #: **M11157**

Product Code: **DSTER-WFI-B-70**

Description: **STER-AHOL 70% WFI Gallon**

Irradiation Batch No.: **2505535**

PRODUCT LOT

Lot # **19-DSTER-910874** EXP **01-AUG-2022**

Lot # EXP

Lot # EXP

Lot # EXP

Specified Dose (K Gy) Min: **25.0** Max: **50.0**

Delivered Dose (K Gy) Min: **27.5** Max: **44.4**

Quantity Irradiated **106 CASES**

A handwritten signature in cursive script, appearing to read "mmfang", is written over a horizontal line.

25Sep2019

Authorized Signature-Quality Assurance

Date

Veltek Associates, Inc.

15 Lee Blvd • Malvern, PA 19355-1234

Tel: (610) 644-8335 Fax: (610) 644-8336

CERTIFICATE OF ANALYSIS



Veltek Associates, Inc.

STER-AHOL70% WFI

LABORATORIES

PRODUCT ANALYSIS DETAILS

Date of Manufacture: 12-AUG-2019

Cat No: **DSTER-WFI-B-70**

Product: **STER-AHOL 70%WFI**

Lot #: **19-DSTER-910874**

Fill Size: **Gallon**

INGREDIENTS

Denatured Ethyl Alcohol
USP Water For Injection

TEST	SPECIFICATIONS	RESULTS
Appearance	Clear colorless	PASS
Assay	68.0-72.0	70.0%
Nonvolatile Residue	<5.0mg	0.0 mg
Specific gravity @ 20 degrees C	0.882-0.892	0.877

m Lang

25Sep2019

Authorized Signature-Quality Assurance

Date

Veltek Associates, Inc.

15 Lee Blvd • Malvern, PA 19355-1234

Tel: (610) 644-8335 Fax: (610) 644-8336

