

MLA Adjustable Volume Models	VistaLab Catalog #	Model	Color	Volume Range	Dispensing Increments	Accuracy*	Precision*
	1132	D-Tipper 50/100/200 μ L selectable pipette	Silver	50–200 μ L		\pm 2.5% at 50 μ L \pm 1.5% at 100 μ L \pm 1.3% at 200 μ L	0.80% at 50 μ L 0.54% at 100 μ L 0.30% at 200 μ L
	1061	Macro 1–5mL adjustable pipette	Blue	1-5mL		\pm 1.0% at 5mL \pm 1.0% at 3mL \pm 1.0% at 1mL	0.60% at 5mL 0.50% at 3mL 0.50% at 1mL
	1134	Macro 2,3,4,5 mL selectable pipette	Blue	2-5mL		\pm 1.0%	0.50%
	1136	Digital pipette 10–50 μ L	Green	10–50 μ L	0.5 μ L	\pm 5% at 10 μ L \pm 3% at 30 μ L \pm 2% at 50 μ L	2.6% at 10 μ L 0.93% at 30 μ L 0.80% at 50 μ L
	1137	Digital pipette 50–200 μ L	Red	50–200 μ L	1.0 μ L	\pm 2% at 50 μ L \pm 1% at 125 μ L \pm 1% at 200 μ L	0.80% at 50 μ L 0.43% at 125 μ L 0.50% at 200 μ L
	1138	Digital pipette 200–1000 μ L	Gray	200–1000 μ L	5.0 μ L	\pm 1% at 200 μ L \pm 1% at 600 μ L \pm 1% at 1000 μ L	0.50% at 200 μ L 1.00% at 600 μ L 0.60% at 1000 μ L

