


The following document has now been written and approved. Information relating to the document supply format is stated below. This document is intended to be circulated within ELGA ONLY.

SECTION 1		FRONT COVER	
Part Number	MANU41196	Version: 01	
Title	PQDIXM1 – PURELAB Quest		
Manual type	Quick Start Manual	<input type="checkbox"/>	Operator Manual <input checked="" type="checkbox"/>





SECTION 2 EXTERNAL CARTON LABELING





SECTION 2-1 Generic position of carton labels & Carton label details

Carton Product Label

Carton UK EPA
Establishment
Number Label




**EPA Est. 091083-GBR-001
PURELAB**


ELGA		
Model	PQDIXXM1	
		
Date of Manufacture	20/09/2019	
Serial No.	PQ000XXXXX	
		
Maximum Pressure		
Inlet	6.0/90	(bar / psi)
Working	6.0/90	(bar / psi)
Electrical Supply	24Vdc-120VA	
	REFER TO OPERATING INSTRUCTIONS SE REFERER AU MANUEL D'OPERATION IN DER BETRIEBSANLEITUNG NACHSEHEN	
ELGA® is the global laboratory water brand of Veolia Water Solutions & Technologies. Tel: +44 203 567 7300, Web: www.elgalabwater.com MADE IN UK		



SECTION 3 PRODUCT LABEL DETAIL


SECTION 3-1 Generic position of units labels & Unit label details.



EPA Est. No. 091083- GBR-001

ELGA 

Model PQDIXM1

 Date of Manufacture 20/09/2019
 Serial No. PQ000XXXX

 Maximum Pressure
 Inlet 6.0/90 (bar / psi)
 Working 6.0/90 (bar / psi)
 Electrical Supply 24Vdc-120VA

 REFER TO OPERATING INSTRUCTIONS
 SE REFERER AU MANUEL D'OPERATION
 IN DER BEDIENUNGSANLEITUNG NACHSEHEN

ELGA® is the global laboratory water brand of
 Veolia Water Solutions & Technologies.
 Tel: +44 203 567 7300, Web: www.elgalabwater.com
 MADE IN UK

AUTHORISATION - RESEARCH AND DEVELOPMENT DEPARTMENT

This document is released on behalf of ELGA and remains the property of VWS (UK) LTD.

Compiled By: <i>Joel Davies</i>	Date: <i>26th September 2019</i>
Approved: <i>[Signature]</i>	Release Date: <i>30th SEPT 19</i>