



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1168-1LP, BDH1168-4LP, BDH1168-5GL, BDH1168-204L, CABDH1168-1LP, CABDH1168-4LP, CABDH1168-5GL, CABDH1168-204L

Catalog Description: Water

Grade: ASTM

Date Current Specification Set: November 24 2009

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Meets ASTM Type II Specifications

Characteristics	Specifications
Electrical Conductivity at 25°C:	Less than 1.0µS/cm
Electrical Resistivity at 25°C:	Greater than or equal to 1.0 megohm-cm
Total Organic Carbon:	Less than 50 µg/L
Sodium (Na):	Less than 5 µg/L
Chlorides (Cl):	Less than 5 µg/L
Total Silica:	Less than 3 µg/L

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor