# MATERIAL SAFETY DATA SHEET

**IDENTITY:** DRIERITE **DATE PREPARED:** 2-21-2011

**DESCRIPTION:** 1/16" TO 1/4" WHITE GRANULES

#### SECTION I

MANUFACTURERS NAME: W.A. HAMMOND DRIERITE CO., LTD. ADDRESS: P.O. BOX 460, 138 DAYTON AVE., XENIA, OHIO 45385

EMERGENCY PHONE NUMBER: 937-376-2927 INFORMATION PHONE NUMBER: 937-376-2927

www.drierite.com

# SECTION II INGREDIENTS

CHEMICAL IDENTITY: % OSHA PEL ACGIH TLV UNITS C.A.S.# TSCA# CALCIUM SULFATE, >99% 15 10 mg//M3 7778-18-9 A739-8876

**ANHYDROUS** 

HAZARDOUS MATERIAL IDENTIFICATION SYSTEMS (HMIS)

**HEALTH**0

FLAMMABILITY
0

REACTIVITY PROTECTIVE EQUIPMENT
1

E

### **SECTION III**

# PHYSICAL/CHEMICAL CHARACTERISTICS

SPECIFIC GRAVITY: (H20-1): 1.87

**SOLUBILITY IN WATER:** 0.25 GRAMS PER LITER

**MELTING POINT: 1450 C DECOMPOSES** 

APPEARANCE: WHITE POWDER OR GRANULES; NO ODOR

### **SECTION IV**

# FIRE AND EXPLOSION HAZARD DATA

**FLASH POINT: NONE** 

**EXTINGUISHER MEDIA:** NOT COMBUSTABLE **SPECIAL FIRE FIGHTING PROCEDURES:** NONE **UNUSUAL FIRE AND EXPLOSION HAZARDS:** NONE

# SECTION V

## **REACTIVITY DATA**

**STABILITY: STABLE** 

INCOMPATIBLE (MATERIALS TO AVOID): STRONG ACIDS HAZARDOUS DECOMPOSITION BYPRODUCTS: SO3 @ 1450 C

HAZARDOUS POLYMERIZATION: WILL NOT OCCUR

#### SECTION VI

# **HEALTH HAZARD DATA**

**EYES: PARTICLES MAY CAUSE IRRITATION** 

SKIN: THIS MATERIAL IS NOT TOXIC. MAY DRY OR IRRITATE SKIN

INHALATION: PERSONS SUBJECTED TO LARGE AMOUNTS OF THIS DUST WILL BE FORCED TO LEAVE AREA BECAUSE OF NUISANCE CONDITIONS SUCH AS COUGHING, SNEEZING,

AND NASAL IRRITATION FROM DUST.

**INGESTION: NOT KNOWN** 

SIGNS AND SYMPTOMS OF OVER EXPOSURE: EYES, NOSE, THROAT, OR RESPIRATORY

**IRRITATION** 

**CARCINOGENICITY OF INGREDIENTS:** 

MATERIAL IARC NTP OSHA

ALL NOT LISTED NOT LISTED NOT LISTED

MEDICAL CONDITIONS GENERALLY AGGRAVATED BY EXPOSURE:

PRE-EXISTING UPPER RESPIRATORY AND LUNG DISEASE SUCH AS, BUT NOT LIMITED TO BRONCHITIS, EMPHYSEMA & ASTHMA

#### **EMERGENCY AND FIRST AID PROCEDURES:**

 $\mbox{\bf EYES:}$  FLUSH WITH WATER. IF IRRITATION CONTINUES OBTAIN MEDICAL ATTENTION

**DUST INHALATION: REMOVE TO FRESH AIR** 

**SKIN:** WASH WITH WATER **INGESTION:** NOT KNOWN

### **SECTION VII**

### SPILL OR LEAKAGE PROCEDURES

# STEPS TO BE TAKEN IN CASE MATERIAL IS RELEASED OR SPILLED:

SWEEP OR VACUUM MATERIAL INTO APPROPRIATE WASTE CONTAINER FOR DISPOSAL. AVOID DUSTING CONDITIONS

WASTE DISPOSAL METHOD: THIS MATERIAL CAN BE DISPOSED OF IN ACCORDANCE WITH PROCEDURES ACCEPTABLE UNDER FEDERAL, STATE, AND LOCAL REGULATIONS. PRECAUTIONS TO BE TAKEN IN HANDLING AND STORING: KEEP CONTAINER CLOSED. STORE IN COOL DRY PLACE. AVOID GENERATING DUST.

# **SECTION VIII**

# **CONTROL MEASURES**

RESPIRATORY PROTECTION: MASK NIOSH/OSHA APPROVED FOR DUST

**VENTILATION:** TO MEET TLV REQUIREMENTS

EYES: SAFETY GLASSES OR GOGGLES

OTHER PROTECTIVE EQUIPMENT: GLOVES OR PROTECTIVE CLOTHING NOT USUALLY

NECESSARY BUT MAY BE DESIRABLE IN SPECIFIC WORK SITUATIONS.

### **SECTION IX**

### **REFERENCES**

U.S. DEPARTMENT OF LABOR - OSHA FORM APPROVED OMB NO. 1218-0072 OSHA HAZARD COMMUNICATION STANDARD-29 CFR 1910. 1200 U.S. GYPSUM CO.

ALTHOUGH THE INFORMATION AND RECOMMENDATION SET FORTH HEREIN ARE PRESENTED IN GOOD FAITH AND BELIEVED TO BE CORRECT AS OF THE DATE HEREOF, THE W.A. HAMMOND DRIERITE CO.,LTD. MAKES NO REPRESENTATION AS TO THE COMPLETENESS OR ACCURACY THEREOF. INFORMATION IS SUPPLIED UPON THE CONDITION THAT THE PERSON RECEIVING IT WILL MAKE THEIR OWN DETERMINATION AS TO ITS SUITABILITY FOR THEIR PURPOSE PRIOR TO USE. IN NO EVENT WILL THE W.A. HAMMOND DRIERITE CO., LTD. BE RESPONSIBLE FOR DAMAGES OF ANY NATURE WHATSOEVER RESULTING FROM THE USE OF OR RELIANCE UPON INFORMATION HEREIN SUPPLIED. NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESSED OR IMPLIED, OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OF OR ANY OTHER NATURE MADE HEREUNDER WITH RESPECT TO INFORMATION OR THE PRODUCT TO WHICH INFORMATION REFERS.