

Steroid Hormone Assay Kits

Researchers have trusted Cayman's collection of steroid hormone ELISAs for decades to study effects in metabolism, inflammation, immune function, osmoregulation, sexual development, and stress responses. Designed with competitive assay principles, each kit offers highly sensitive detection in the low pg/ml range, reliable and reproducible results, and flexible use with multiple sample types.

Androgens

Product Name	Cayman Item No.	VWR Cat. No.
Testosterone ELISA Kit	582701	75816-942 75816-940
11-keto Testosterone ELISA Kit	582751	75816-946 75816-944

Oestrogens

Product Name	Cayman Item No.	VWR Cat. No.
Estradiol ELISA Kit	582251	75817-380 75817-378
Estriol ELISA Kit	582281	75816-934 75816-932
2-Methoxyestradiol ELISA Kit	582261	75817-384 75817-382

Progestins

Product Name	Cayman Item No.	VWR Cat. No.
Progesterone ELISA Kit	582601	75816-938 75816-936

Glucocorticoids


Product Name	Cayman Item No.	VWR Cat. No.
Corticosterone ELISA Kit	501320	75816-864 75816-862
Cortisol ELISA Kit	500360	75816-780 75816-778
Cortisol Express ELISA Kit	500370	75816-784 75816-782

Vitamin D

Product Name	Cayman Item No.	VWR Cat. No.
Vitamin D ELISA Kit	501050	75816-842 75816-840

Why Choose Cayman's Steroid Assay Kits?

Cayman's steroid assays undergo rigorous quality testing to deliver the sensitivity and specificity needed to detect biologically significant analyte levels. These kits were developed to give the best sensitivity possible, while maintaining good signal-to-noise ratio, and report sensitivity as the 80% B/B₀ point on the standard curve. This point of quantitation can be trusted to give you reliable, reproducible data.



For additional product or technical information, visit www.vwr.com, call 1.800.932.5000, or contact your VWR representative.