

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

**SECTION 1. IDENTIFICATION**

Product name : Reagent Alcohol

Number : 000000011364

Product Use Description : Solvent

Manufacturer or supplier's details : Honeywell International Inc.  
1953 South Harvey Street  
Muskegon, MI 49442

For more information call : 1-800-368-0050  
+1-231-726-3171  
  
(Monday-Friday, 9:00am-5:00pm)

**In case of emergency call : Medical: 1-800-498-5701 or +1-303-389-1414**  
: **Transportation (CHEMTREC): 1-800-424-9300 or**  
: **+1-703-527-3887**  
:  
: (24 hours/day, 7 days/week)

**SECTION 2. HAZARDS IDENTIFICATION****Emergency Overview**

Form : liquid, clear

Color : colourless

Odor : mild alcohol-like

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

**Classification of the substance or mixture**

Classification of the substance or mixture : Flammable liquids, Category 2  
Eye irritation, Category 2A  
Specific target organ toxicity - single exposure, Category 1,  
Eyes, Nervous system, Systemic toxicity

**GHS Label elements, including precautionary statements**

Symbol(s)



Signal word

: Danger

Hazard statements

: Highly flammable liquid and vapour.  
Causes serious eye irritation.  
Causes damage to organs.

Precautionary statements

**Prevention:**  
Keep away from heat/sparks/open flames/hot surfaces. No smoking.  
Keep container tightly closed.  
Ground/bond container and receiving equipment.  
Use explosion-proof electrical/ ventilating/ lighting/ equipment.  
Use only non-sparking tools.  
Take precautionary measures against static discharge.  
Do not breathe dust/ fume/ gas/ mist/ vapours/ spray.  
Wash skin thoroughly after handling.  
Do not eat, drink or smoke when using this product.  
Wear protective gloves/ protective clothing/ eye protection/ face protection.

**Response:**  
IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower.  
IF IN EYES: Rinse cautiously with water for several minutes.

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

Remove contact lenses, if present and easy to do. Continue rinsing.  
IF exposed: Call a POISON CENTER or doctor/ physician.  
If eye irritation persists: Get medical advice/ attention.  
In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction.

**Storage:**

Store in a well-ventilated place. Keep cool.  
Store locked up.

**Disposal:**

Dispose of contents/ container to an approved waste disposal plant.

**Carcinogenicity**

ACGIH: Ethanol 64-17-5  
A3: Confirmed animal carcinogen

**SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS**

Chemical nature : Mixture

Chemical name	CAS-No.	Concentration
Ethanol	64-17-5	90.00 %
Methanol	67-56-1	5.00 %
Isopropanol	67-63-0	5.00 %

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

**SECTION 4. FIRST AID MEASURES**

- Inhalation : Call a physician immediately. Remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Use oxygen as required, provided a qualified operator is present.
- Skin contact : Wash off immediately with plenty of water for at least 15 minutes. Take off contaminated clothing and shoes immediately. Wash contaminated clothing before re-use. Call a physician.
- Eye contact : Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Call a physician.
- Ingestion : Call a physician immediately. Do NOT induce vomiting. Immediate medical attention is required. Never give anything by mouth to an unconscious person.

**Notes to physician**

- Indication of immediate medical attention and special treatment needed, if necessary : Treat symptomatically.

**SECTION 5. FIREFIGHTING MEASURES**

- Suitable extinguishing media : Alcohol-resistant foam  
Carbon dioxide (CO<sub>2</sub>)  
Dry chemical  
Cool closed containers exposed to fire with water spray.
- Unsuitable extinguishing media : Do not use a solid water stream as it may scatter and spread fire.
- Specific hazards during firefighting : Flammable.  
Vapours may form explosive mixtures with air.  
Vapours are heavier than air and may spread along floors.

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

Vapors may travel to areas away from work site before igniting/flashing back to vapor source.  
In case of fire hazardous decomposition products may be produced such as:  
Carbon monoxide  
Carbon dioxide (CO<sub>2</sub>)  
Formaldehyde

Special protective equipment for firefighters : Wear self-contained breathing apparatus and protective suit.

**SECTION 6. ACCIDENTAL RELEASE MEASURES**

- Personal precautions, protective equipment and emergency procedures : Wear personal protective equipment.  
Immediately evacuate personnel to safe areas.  
Keep people away from and upwind of spill/leak.  
Ensure adequate ventilation.  
Remove all sources of ignition.  
Do not swallow.  
Do not breathe vapours or spray mist.  
Avoid contact with skin, eyes and clothing.
- Environmental precautions : Prevent further leakage or spillage if safe to do so.  
Prevent product from entering drains.  
Discharge into the environment must be avoided.  
Do not flush into surface water or sanitary sewer system.  
Do not allow run-off from fire fighting to enter drains or water courses.
- Methods and materials for containment and cleaning up : Ventilate the area.  
No sparking tools should be used.  
Use explosion-proof equipment.  
Contain spillage, and then collect with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see section 13).

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

**SECTION 7. HANDLING AND STORAGE****Handling**

Precautions for safe handling : Wear personal protective equipment.  
Use only in well-ventilated areas.  
Keep container tightly closed.  
Do not smoke.  
Do not swallow.  
Do not breathe vapours or spray mist.  
Avoid contact with skin, eyes and clothing.

Advice on protection against fire and explosion : Keep away from fire, sparks and heated surfaces.  
Take precautionary measures against static discharges.  
Ensure all equipment is electrically grounded before beginning transfer operations.  
Use explosion-proof equipment.  
Keep product and empty container away from heat and sources of ignition.  
No sparking tools should be used.  
No smoking.

**Storage**

Conditions for safe storage, including any incompatibilities : Store in area designed for storage of flammable liquids. Protect from physical damage.  
Keep containers tightly closed in a dry, cool and well-ventilated place.  
Containers which are opened must be carefully resealed and kept upright to prevent leakage.  
Keep away from heat and sources of ignition.  
Keep away from direct sunlight.  
Store away from incompatible substances.  
Container hazardous when empty.  
Do not pressurize, cut, weld, braze, solder, drill, grind or expose containers to heat or sources of ignition.

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

**SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION**

- Protective measures : Ensure that eyewash stations and safety showers are close to the workstation location.
- Engineering measures : Use with local exhaust ventilation.  
Prevent vapour buildup by providing adequate ventilation during and after use.
- Eye protection : Do not wear contact lenses.  
Wear as appropriate:  
Safety glasses with side-shields  
If splashes are likely to occur, wear:  
Goggles or face shield, giving complete protection to eyes
- Hand protection : Solvent-resistant gloves  
Gloves must be inspected prior to use.  
Replace when worn.
- Skin and body protection : Wear as appropriate:  
Solvent-resistant apron  
Flame retardant antistatic protective clothing.  
If splashes are likely to occur, wear:  
Protective suit
- Respiratory protection : In case of insufficient ventilation wear suitable respiratory equipment.  
For rescue and maintenance work in storage tanks use self-contained breathing apparatus.  
Use NIOSH approved respiratory protection.
- Hygiene measures : When using, do not eat, drink or smoke.  
Wash hands before breaks and immediately after handling the product.  
Keep working clothes separately.  
Remove and wash contaminated clothing before re-use.  
Do not swallow.  
Do not breathe vapours or spray mist.  
Avoid contact with skin, eyes and clothing.

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

## Exposure Guidelines

Components	CAS-No.	Value	Control parameters	Update	Basis
Ethanol	64-17-5	STEL : Short term exposure limit	(1,000 ppm)	2009	ACGIH:US. ACGIH Threshold Limit Values
Ethanol	64-17-5	REL : Recommended exposure limit (REL):	1,900 mg/m3 (1,000 ppm)	2005	NIOSH/GUIDE:US. NIOSH: Pocket Guide to Chemical Hazards
Ethanol	64-17-5	PEL : Permissible exposure limit	1,900 mg/m3 (1,000 ppm)	02 2006	OSHA_TRANS:US. OSHA Table Z-1 Limits for Air Contaminants (29 CFR 1910.1000)
Ethanol	64-17-5	TWA : Time weighted average	1,900 mg/m3 (1,000 ppm)	1989	Z1A:US. OSHA Table Z-1-A (29 CFR 1910.1000)
Methanol	67-56-1	SKIN_DES : Skin designation:	Can be absorbed through the skin.	2008	ACGIH:US. ACGIH Threshold Limit Values
Methanol	67-56-1	STEL : Short term exposure limit	(250 ppm)	2008	ACGIH:US. ACGIH Threshold Limit Values



## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

Methanol	67-56-1	TWA : Time weighted average	(200 ppm)	2008	ACGIH:US. ACGIH Threshold Limit Values
Methanol	67-56-1	STEL : Short term exposure limit	325 mg/m3 (250 ppm)	2005	NIOSH/GUIDE:US. NIOSH: Pocket Guide to Chemical Hazards
Methanol	67-56-1	SKIN_DE S : Skin designati on:	Can be absorbed through the skin.	2005	NIOSH/GUIDE:US. NIOSH: Pocket Guide to Chemical Hazards
Methanol	67-56-1	REL : Recomm ended exposure limit (REL):	260 mg/m3 (200 ppm)	2005	NIOSH/GUIDE:US. NIOSH: Pocket Guide to Chemical Hazards
Methanol	67-56-1	PEL : Permissi ble exposure limit	260 mg/m3 (200 ppm)	02 2006	OSHA_TRANS:US. OSHA Table Z-1 Limits for Air Contaminants (29 CFR 1910.1000)
Methanol	67-56-1	SKIN_FI NAL : Skin designati on (Final Rule Limit applies):	Can be absorbed through the skin.	1989	Z1A:US. OSHA Table Z-1-A (29 CFR 1910.1000)

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

Methanol	67-56-1	STEL : Short term exposure limit	325 mg/m <sup>3</sup> (250 ppm)	1989	Z1A:US. OSHA Table Z-1-A (29 CFR 1910.1000)
Methanol	67-56-1	TWA : Time weighted average	260 mg/m <sup>3</sup> (200 ppm)	1989	Z1A:US. OSHA Table Z-1-A (29 CFR 1910.1000)
Isopropanol	67-63-0	TWA : Time weighted average	(200 ppm)	2008	ACGIH:US. ACGIH Threshold Limit Values
Isopropanol	67-63-0	STEL : Short term exposure limit	(400 ppm)	2008	ACGIH:US. ACGIH Threshold Limit Values
Isopropanol	67-63-0	REL : Recomm ended exposure limit (REL):	980 mg/m <sup>3</sup> (400 ppm)	2005	NIOSH/GUIDE:US. NIOSH: Pocket Guide to Chemical Hazards
Isopropanol	67-63-0	STEL : Short term exposure limit	1,225 mg/m <sup>3</sup> (500 ppm)	2005	NIOSH/GUIDE:US. NIOSH: Pocket Guide to Chemical Hazards
Isopropanol	67-63-0	PEL : Permissi ble exposure limit	980 mg/m <sup>3</sup> (400 ppm)	02 2006	OSHA_TRANS:US. OSHA Table Z-1 Limits for Air Contaminants (29 CFR 1910.1000)

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

Isopropanol	67-63-0	TWA : Time weighted average	980 mg/m3 (400 ppm)	1989	Z1A:US. OSHA Table Z-1-A (29 CFR 1910.1000)
-------------	---------	--------------------------------------	------------------------	------	---------------------------------------------------

Isopropanol	67-63-0	STEL : Short term exposure limit	1,225 mg/m3 (500 ppm)	1989	Z1A:US. OSHA Table Z-1-A (29 CFR 1910.1000)
-------------	---------	----------------------------------------------	--------------------------	------	---------------------------------------------------

## SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Physical state	: liquid, clear
Color	: colourless
Odor	: mild alcohol-like
Odor threshold	: Note: no data available
pH	: Note: Not applicable
Melting point/range	: -114.1 °C Note: The physical data is that of the main component.
Boiling point/boiling range	: 78.32 °C Note: The physical data is that of the main component.
Flash point	: 59 °F (15 °C) Method: closed cup Note: The physical data is that of the main component.
Evaporation rate	: ca. 3 Method: Compared to Butyl acetate.

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

Lower explosion limit	: 3.3 %(V) Note: The physical data is that of the main component.
Upper explosion limit	: 19 %(V) Note: The physical data is that of the main component.
Vapor pressure	: 59.5 hPa at 20 °C(68 °F)Note: The physical data is that of the main component.
Vapor density	: 1.6 Note: (Air = 1.0), The physical data is that of the main component.
Density	: 0.7890 g/cm <sup>3</sup> at 20 °C 0.7847 g/cm <sup>3</sup> at 25 °C
Water solubility	: Note: completely soluble
Partition coefficient: n-octanol/water	: Note: no data available
Ignition temperature	: Note: not determined
Decomposition temperature	: Note: no data available
Viscosity, dynamic	: Note: no data available
Viscosity, kinematic	: Note: no data available

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

**SECTION 10. STABILITY AND REACTIVITY**

Chemical stability	: Stable under recommended storage conditions.
Possibility of hazardous reactions	: Hazardous polymerisation does not occur.
Conditions to avoid	: Heat, flames and sparks. Keep away from direct sunlight.
Incompatible materials	: Strong oxidizing agents Potassium superoxide Bromine Pentafluoride Acetyl bromide Acetyl chloride Platinum Sodium
Hazardous decomposition products	: In case of fire hazardous decomposition products may be produced such as: Carbon monoxide Carbon dioxide (CO <sub>2</sub> ) Formaldehyde

**SECTION 11. TOXICOLOGICAL INFORMATION**

Acute oral toxicity Ethanol	: LD50: 7,060 mg/kg Species: Rat
Methanol	: LD50: 5,628 mg/kg Species: Rat

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

## Acute inhalation toxicity

Ethanol : LC50: 20000 ppm  
Exposure time: 10 h  
Species: Rat

Methanol : LC50: 64000 ppm  
Exposure time: 4 h  
Species: Rat

## Acute dermal toxicity

Methanol : LD50: 15,800 mg/kg  
Species: Rabbit

## Skin irritation

Ethanol : Species: Rabbit  
Result: Irritating to skin.  
Exposure time: 24 h

Methanol : Species: Rabbit  
Classification: irritating  
Exposure time: 24 h

## Eye irritation

Ethanol : Species: Rabbit  
Result: Irritating to eyes.  
Exposure time: 24 h

Methanol : Species: rabbit eye  
Classification: irritating

## Repeated dose toxicity

Methanol : Species: Rat  
Application Route: Inhalation  
Test substance: Methanol  
Developmental Toxicity  
NOAEL (maternal toxicity)  
10,000 ppm  
NOAEL (developmental toxicity)

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

5,000 ppm  
Skeletal and visceral malformations.

Genotoxicity in vitro  
Methanol : Note: In vitro tests did not show mutagenic effects

Genotoxicity in vivo  
Methanol : Note: In vivo tests did not show mutagenic effects

Further information : Note: Confirmed animal carcinogen with unknown relevance to humans.

**SECTION 12. ECOLOGICAL INFORMATION**

Toxicity to fish  
Ethanol : LC0: 8,140 mg/l  
Exposure time: 48 h  
Species: Leuciscus idus (Golden orfe)

flow-through test  
LC50: 12,900 mg/l  
Exposure time: 96 h  
Species: Oncorhynchus mykiss (rainbow trout)

LC50: 14,200 mg/l  
Exposure time: 96 h  
Species: Pimephales promelas (fathead minnow)

Methanol : LC50: 29,400 mg/l  
Exposure time: 96 h  
Species: Pimephales promelas (fathead minnow)

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

## Toxicity to daphnia and other aquatic invertebrates

Ethanol : EC50: 9,268 mg/l  
Exposure time: 48 h  
Species: Daphnia magna (Water flea)

EC50: 10,800 mg/l  
Exposure time: 24 h  
Species: Daphnia magna (Water flea)

Methanol : LC50: 10,000 mg/l  
Exposure time: 24 h  
Species: Daphnia (water flea)

## Toxicity to algae

Ethanol : LC0: 5,000 mg/l  
Species: Scenedesmus quadricauda (Green algae)

## Toxicity to bacteria

Ethanol : LC0: 6,500 mg/l  
Species: Pseudomonas putida

EC50: 35,470 mg/l  
Exposure time: 5 min  
Species: Photobacterium phosphoreum

EC50: 34,634 mg/l  
Exposure time: 30 min  
Species: Photobacterium phosphoreum

Methanol : EC50: 43,000 mg/l  
Exposure time: 5 min  
Species: Photobacterium phosphoreum

EC50: 40,000 mg/l  
Exposure time: 15 min  
Species: Photobacterium phosphoreum

EC50: 39,000 mg/l  
Exposure time: 25 min  
Species: Photobacterium phosphoreum



**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

Biodegradability  
Isopropanol : Biochemical Oxygen Demand (BOD) Biochemical oxygen demand within 5 days  
Value: 58 %

**Further information on ecology****SECTION 13. DISPOSAL CONSIDERATIONS**

Disposal methods : Observe all Federal, State, and Local Environmental regulations.

**SECTION 14. TRANSPORT INFORMATION**

**DOT** UN/ID No. : UN 1987  
Proper shipping name : ALCOHOLS, N.O.S.  
(Ethanol, Methanol, Isopropanol)  
Class : 3  
Packing group : II  
Hazard Labels : 3

**IATA** UN/ID No. : UN 1987  
Description of the goods : ALCOHOLS, N.O.S.  
(Ethanol, Methanol, Isopropanol)  
Class : 3  
Packaging group : II  
Hazard Labels : 3  
Packing instruction (cargo aircraft) : 364  
Packing instruction (passenger aircraft) : 353  
Packing instruction (passenger aircraft) : Y341

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

<b>IMDG</b>	UN/ID No.	: UN 1987
	Description of the goods	: ALCOHOLS, N.O.S. (ETHANOL, METHANOL, ISOPROPANOL)
	Class	: 3
	Packaging group	: II
	Hazard Labels	: 3
	EmS Number	: F-E, S-D
	Marine pollutant	: no

**SECTION 15. REGULATORY INFORMATION****Inventories**

US. Toxic Substances Control Act : On TSCA Inventory

Australia. Industrial Chemical (Notification and Assessment) Act : On the inventory, or in compliance with the inventory

Canada. Canadian Environmental Protection Act (CEPA). Domestic Substances List (DSL) : All components of this product are on the Canadian DSL

Japan. Kashin-Hou Law List : On the inventory, or in compliance with the inventory

Korea. Existing Chemicals Inventory (KECI) : On the inventory, or in compliance with the inventory

Philippines. The Toxic Substances and Hazardous and Nuclear Waste Control Act : On the inventory, or in compliance with the inventory

China. Inventory of Existing Chemical Substances : On the inventory, or in compliance with the inventory

New Zealand. Inventory of Chemicals (NZIoC), as

## Reagent Alcohol

## LP090-4P

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

published by ERMA New  
Zealand

**National regulatory information**

US. EPA CERCLA Hazardous Substances (40 CFR 302) : The following component(s) of this product is/are subject to release reporting under 40 CFR 302 when release exceeds the Reportable Quantity (RQ):

Reportable quantity: 5000 lbs  
: Methanol 67-56-1

**SARA 302 Components** : No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

**SARA 313 Components** : The following components are subject to reporting levels established by SARA Title III, Section 313:

: Methanol 67-56-1  
: Isopropanol 67-63-0

**SARA 311/312 Hazards** : Fire Hazard  
Acute Health Hazard  
Chronic Health Hazard

**California Prop. 65** :



**WARNING:** This product can expose you to chemicals, listed below, known to the State of California to cause birth defects or other reproductive harm. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Methanol 67-56-1

**Massachusetts RTK** : Ethanol 64-17-5  
: Methanol 67-56-1  
: Isopropanol 67-63-0

**New Jersey RTK** : Ethanol 64-17-5  
: Methanol 67-56-1  
: Isopropanol 67-63-0

**Reagent Alcohol****LP090-4P**

Version 2.1

Revision Date 08/16/2018

Print Date 10/17/2019

<b>Pennsylvania RTK</b>	: Ethanol	64-17-5
	: Methanol	67-56-1
	: Isopropanol	67-63-0

**SECTION 16. OTHER INFORMATION**

	<b>HMIS III</b>	<b>NFPA</b>
Health hazard	: 2*	0
Flammability	: 3	3
Physical Hazard	: 0	
Instability	:	0

\* - Chronic health hazard

Hazard rating and rating systems (e.g. HMIS® III, NFPA): This information is intended solely for the use of individuals trained in the particular system.

**Further information**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text. Final determination of suitability of any material is the sole responsibility of the user. This information should not constitute a guarantee for any specific product properties.

Changes since the last version are highlighted in the margin. This version replaces all previous versions.

Previous Issue Date: 03/12/2018

Prepared by Honeywell Performance Materials and Technologies Product Stewardship Group