

Product Name:	HCUCFS0120
Description:	<p>VWR Catalog No: 10819-662 VWR Countertop Laboratory Freezer Freestanding 1.7 Cu. Ft. Capacity One solid swing door, right hinged* Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in glycol bottle Adjustable operating temperature range: -15°C – -25°C 2 adjustable shelves with guard rail to prevent product from being pushed against back wall Bottom footed shelf to allow air circulation Keyed door lock Magnetic door gasket for positive seal Leveling legs Integrated door handle Probe access port (3/8") Powder coat exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Exterior dimensions: 23 3/4 "W x 24" D x 20 3/8"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Manual defrost 115V, 60 Hz, 1 Amp, 1/8 HP UL/C-UL listed Shipping weight: 106 lbs. *This product is available left hinged, product number HCUCBI0404LH</p>