



Material	75164-300
Material description	Potassium sodium tartrate
Grade	ACS
CAS Number	6381-59-5
Molecular formula	KOOC(CHOH)2COONa·4H2O
Molecular mass	282.22

Characteristics

Specifications

Assay	99.0 --> 102.0 %
pH (25°C; 5 %)	6.0 --> 8.5
Heavy metals (as Pb)	Max. 5 ppm
Insoluble matter	Max. 0.005 %
Cl (Chloride)	Max. 0.001 %
NH4 (Ammonium)	Max. 0.002 %
PO4 (Phosphate)	Max. 0.002 %
SO4 (Sulphate)	Max. 0.005 %
Ca (Calcium)	Max. 0.005 %
Fe (Iron)	Max. 0.001 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.