



VWR ANALYTICAL

## Product Specification

**Catalog Number:** BDH1115-1LP, BDH1115-4LP, CABDH1115-1LP, CABDH1115-4LP

**Catalog Description:** Dimethyl Sulfoxide

**Grade:** ACS

**Date Current Specification Set:** April 19 2006

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

### Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance	Clear, colorless liquid
Water (K.F.):	Less than 0.10%
Residue	Less than 0.01% <b>(0.0010g/91mL)</b>
GC Purity:	Greater than 99.9%
Color (APHA):	Less than 10
Titration acid	Less than 1 $\mu$ eq/g (Not more than 5.0mL of of 0.01N sodium hydroxide)

Signed on behalf of VWR:

Neal Fox  
Quality Control Supervisor