

Product Specification

Meets U.S.P Requirements, GMP Manufactured Product

Test	Specification
USP – Assay (Na ₂ SO ₄) (dried basis)	≥ 99.0 %
USP – Acidity or Alkalinity (mL)	≤ 0.50
USP – Chloride (Cl)	≤ 0.005 %
USP – Identification A	Passes Test
USP – Identification B	Passes Test
USP – Insoluble Matter	≤ 0.01 %
USP – Loss on Drying at 105°C	≤ 0.5 %
Thru U.S. No. 50 Sieve, For Informa	%

Bulk Pharmaceutical Chemical
GMP Manufactured Product

CAUTION: For Manufacturing, processing or repackaging

No Class 1,2,3 or other solvents are used or produced in the manufacturing or purification of the product.

Metallic Residues: No metal catalysts or metal reagents, as defined by EMA Guideline EMEA/CHMP/SWP/4446/2000 , are used in the production of this material.

Appearance – Colorless to slightly colored crystalline solid.

Packaging Site: Paris Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality