


| | | | | | |
|---------------------|---|---------------|---------------|---------------|---|
| PDS No. 78190x | PRODUCT DATA SHEET | | | | Page 1 of 1 |
| Revision 03 | 384 Well Microplate, PS, F-Bottom, Non-Binding | | | |  |
| | Greiner Item-No. 78190x | | | | |
| Valid for Item-No.: | 781900 | 781901 | 781903 | 781904 | 781906 |

| 1. | Description / Specification | |
|-----|-----------------------------|---|
| 1.1 | Description | PS Microplate, 384 well, rounded square well design, alphanumeric well coding, protein-repellent Non-Binding-Treatment 781900: solid F-bottom (flat) 781901: solid F-bottom (flat) 781903: µClear® F-bottom (flat) 781904: solid F-bottom (flat) 781906: µClear® F-bottom (flat) |
| 1.2 | Dimensions | Length: 127,76 +/- 0,15 mm Width: 85,48 +/- 0,15 mm Curvature: ≤ 150 µm |
| 1.3 | Volume per well | Total volume: 131 µl Working volume: 15 -110 µl |
| 1.4 | Material / Resin | Modified PS (Polystyrene), free of heavy metal |
| 1.5 | Colour | 781900: black 781901: clear 781903: white plate, clear bottom 781904: white 781906: black plate, clear bottom |
| 1.6 | Sterilization | No |
| 1.7 | Quality Control | - <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Other Information | For single use only |

| 2. | Features | |
|-----|--------------------------|---|
| 2.1 | Basic features | Free of detectable DNase/RNase, human DNA and pyrogens |
| 2.2 | Temperature range | -20°C to +60°C |
| 2.3 | Autoclavability | No |
| 2.4 | Centrifugation, max. RCF | 4800 x g: swinging-bucket rotor |
| 2.5 | Chemical Resistance | See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins |
| 2.6 | Shelf life | 4 years after month of production |
| 2.7 | Other Information | - |

| 3. | Packaging | |
|-----|-------------------|--|
| 3.1 | Pieces / Bag | 10 |
| 3.2 | Pieces / Box | 40 |
| 3.3 | Lot-No. | E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) |
| 3.4 | Other Information | Certificate of Quality |

| 4. | Other Information |
|----|-------------------|
| | - |

Data Sheet subject to change without notice!

| Prior Issue | Drawn | Approved | Released | CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|----------------------|-----------------------|----------------------|--|
| Revision 02 | Date 24 June 2013 | Date 24 June 2013 | Date 24 June 2013 | |
| Date 10.12.2009 | Name S. Kaelberer | Name Dr. R. Heller | Name A. Schulz | |