

Product Name:	VWR-HC-SSP-23FA
Description:	<p>VWR Catalog Number: 76514-214  VWR Plus Series Stainless Steel Freezer  23 Cu. Ft. capacity  One swing solid door, self-closing, right hinged, not reversible  3 ea. Adjustable Shelves  Digital display microprocessor temperature controller  Adjustable operating temperature: -15°C to -25°C  Audible and visual Hi/Lo temperature alarm  2 probes; 1 in air and 1 in sample bottle  Probe attached to evaporator for defrost termination sensing  Sensor error alarm  Continuous product (simulated) temperature display  Remote alarm contacts  Forced draft air circulation  Keyed door lock  Magnetic door gaskets for positive seal  Integrated door handle  LED interior lights are safety shielded and switch controlled  Probe access port (3/4")  Stainless steel exterior  High density urethane foam cabinet insulation  White-powder coated aluminum interior with 304 stainless steel floor  Easily cleaned smooth interior  Hydrocarbon natural refrigerant (R290)  2 year parts and labor warranty, plus an additional 3 year  compressor parts warranty  Auto defrost  4 swivel casters (5 3/4")  Exterior dimensions: 26 7/8" W x 35 3/4" D x 83 1/4" H with casters;  77 1/2" H without casters  115V, 60 Hz, 8 Amps, 3/4 HP  ETL listed  Shipping weight: 327 lbs.</p>