

## Product Specification

Test	Specification
Normality	4.975 – 5.025
Sulfate (SO <sub>4</sub> )	≤ 1 ppm
Ammonium (NH <sub>4</sub> )	≤ 3 ppm
Trace Impurities – Heavy Metals (as Pb)	≤ 0.5 ppm
Trace Impurities – Iron (Fe)	≤ 0.1 ppm
Appearance	Passes Test

For Laboratory, Research, or Manufacturing Use  
Standardization at 25°C traceable to NIST Standard Reference Material.

Packaging Site: Paris Mfg Ctr & DC