

|               |  |
|---------------|--|
| Product Name: | GPR231SSS/0  |
| Description:  | <p>Corepoint Scientific™ General Purpose Stainless Steel Refrigerator<br/> 23 Cu. Ft. Capacity<br/> One swing solid door, self-closing, right hinged, not reversible<br/> 3 ea. Adjustable Shelves<br/> Digital display microprocessor temperature controller<br/> Adjustable operating temperature: 1°C to 10°C (34°F to 50°F)<br/> Audible and visual Hi/Lo temperature alarm<br/> 2 probes; 1 in air and 1 in sample bottle<br/> Sensor error alarm<br/> Continuous product (simulated) temperature display<br/> Remote alarm contacts<br/> Forced draft air circulation<br/> Keyed door lock<br/> Magnetic door gaskets for positive seal<br/> Integrated door handle<br/> LED interior lights are safety shielded and switch controlled<br/> Probe access port (3/4")<br/> Stainless steel exterior<br/> White-powder coated aluminum interior with 304 stainless steel floor<br/> Easily cleaned smooth interior<br/> Hydrocarbon natural refrigerant (R290)<br/> 2 year parts and labor warranty, plus an additional 3 year<br/> compressor parts warranty<br/> 4 swivel casters (5 3/4")<br/> Exterior dimensions: 26 7/8" W x 33 3/4" D x 83 1/4" H with casters;<br/> 77 1/2" H without casters<br/> Cycle defrost<br/> 115V, 60 Hz, 3 Amps, 1/5 HP<br/> ETL listed<br/> Shipping weight: 311 lbs.</p> |