



Material	46644-700
Material description	Hydrogen peroxide 30%
Grade	ACS
CAS Number	7722-84-1
Molecular formula	H2O2
Molecular mass	34.01

## Characteristics

## Specifications

Assay	29.0 --> 32.0 %
Color (Pt-Co)	Max. 10 APHA
Titration acid	Max. 0.0006 meq/g
Heavy metals (as Pb)	Max. 1 ppm
Residue after evaporation	Max. 0.002 %
Cl (Chloride)	Max. 2 ppm
NH4 (Ammonium)	Max. 5 ppm
NO3 (Nitrate)	Max. 2 ppm
PO4 (Phosphate)	Max. 2 ppm
SO4 (Sulphate)	Max. 5 ppm
Fe (Iron)	Max. 0.5 ppm

## Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine  
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,  
QC, Canada

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.