

Controlled Environments + Industrial Products





Puritan® CONTROLLED ENVIRONMENTS + INDUSTRIAL

Puritan's single-use products for industrial, controlled and critical environments are manufactured in our certified and CE-compliant facility, providing reliable, consistent performance to ensure the quality of your end product.

Our swabs are designed for use in specific environments such as precision manufacturing, fiber optics, electronics, semiconductors, and biotechnology. For even tighter control, we offer critical environment swabs for the cleanroom.

Our line includes foam, cotton, knitted polyester, ESD (electrostatic discharge), anti-static, and chemical-resistant swabs.

Polyester and microfiber tips, soft but strong

Each swab is made to help meet your QC goals in a specific, controlled environment.

Puritan's polyester and microfiber swabs are easy on delicate components yet strong enough to stand up to industrial use. Manufactured in Puritan's cleanroom, these swabs are protected from contamination. Personnel wear smocks, gloves, hair covers and masks.

SMALL TIP
Ideal for fine assembly tasks.

FLEXIBLE HANDLE
Makes it easy to reach areas that need to be swabbed. Polypropylene shaft is also resistant to breakdown from chemical content.

THERMAL BONDING
Eliminates the need for adhesives that could transfer to the surface being cleaned.

LOW IN NVR
Knitted polyester is ideal for removing residues and particulates.

FEATURED PRODUCTS:

Left: PurSwab® knitted polyester-tipped applicator (#89133-794 / #3135)

Middle: PurSwab® knitted polyester-tipped applicator (#89133-800 / #3655)









Right: PurSwab® microfiber-tipped applicator (#890897-856 / #3677)

Did you know?

Microfibers are half the diameter of fine silk fiber, one-quarter the diameter of fine wool, and **100 times finer than human hair.**



PURSWAB® KNITTED POLYESTER-TIPPED APPLICATORS

Knitted polyester is non-abrasive and lint-free. In addition, the polyester can easily stand up to IPA, acetone and other solvents commonly used in cleanrooms and other critical environments.

	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	89133-794	3135	Small rigid tip, polypropylene handle	11.989 x 4.496 mm 0.472 x 0.177"	68.504 mm 2.697"	10/100, 1 lb
	89133-802	3140	Thin flexible tip, polypropylene handle	18.999 x 7.29 mm 0.748 x 0.287"	162 mm 6.378"	10/100, 4 lb
	82030-704	3600	Large paddle tip, polypropylene handle	28.575 x 9.906 mm 1.125 x 0.39"	157.099 mm 6.185"	20/50, 3 lb
	82030-706	3605	Symmetrical round tip, polypropylene handle	19.837 x 5.156 mm 0.781 x 0.203"	152.4 mm 6"	10/100, 3 lb
	82030-710	3615	Half-round tip, acetal handle	14.605 x 3.175 mm 0.575 x 0.125"	68.275 mm 2.688"	10/100, 2 lb
	82030-714	3625	Symmetrical round tip, acetal handle	10.312 x 3.175 mm 0.406 x 0.125"	71.755 mm 2.825"	10/100, 2 lb
	89133-800	3655	Large flat paddle tip, polypropylene handle	25.908 x 13.97 mm 1.02 x 0.55"	127.787 mm 5.031"	10/100, 5 lb
	89133-806	3676	Flexible paddle tip, polypropylene handle	16.789 x 6.807 mm 0.661 x 0.268"	154.838 mm 6.096"	20/50, 2 lb

PURSWAB® MICROFIBER-TIPPED APPLICATORS

Microfiber tips are soft and non-abrasive, with very low particulates. These swabs are resistant to normal cleaning solvents, which makes them ideal for removing contaminants and excess materials in fine assembly environments.

	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	89133-808	3130	Small flexible paddle tip, polypropylene handle	10.008 x 3.429 mm 0.394 x 0.135"	69.012 mm 2.717"	10/100, 1 lb
	89087-856	3677	Flexible paddle tip, polypropylene handle	16.789 x 6.807 mm 0.661 x 0.268"	154.838 mm 6.096"	10/100, 2 lb

PurSwab® foam swabs

They are soft, non-abrasive, and residue-free.

Choose from configurations that offer rigid or flexible handles – with foam tips designed for a range of standard and specialty purposes.



#10806-032 / #1805-PF RECT

REGULAR FOAM
For general cleaning, suitable for applying solutions and solvents.



#89087-862 / #1621-PF ESD

ESD FOAM
Small point for precision cleaning. ESD plastic handle.



#10806-028 / #1806-PCFL

FOAM OVER COTTON
Cotton core for increased absorbency, foam mitt for lint-sensitive applications.



#10806-154 / #1806-PF AS

ANTI-STATIC FOAM
Abrasive for cleaning surfaces.



#10806-144 / #1806-WCF CR

CHEMICAL RESISTANT
Tip will not swell or break apart with most solvents.

FEATURED PRODUCTS:

A selection of products from our PurSwab® line

PURSWAB® FOAM-TIPPED APPLICATORS

Thermal edge seals and thermal bonding to the handle help ensure the sponge-like foam will clean up unwanted particles without releasing any into your critical environment.







	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	89133-790	1130	Pointed tip, polypropylene handle	11.506 x 3.505 mm 0.453 x 0.138"	69.494 mm 2.736"	10/100, 2 lb
	89133-792	1135	Thin flexible paddle tip, polypropylene handle	14.224 x 3.962 mm 0.56 x 0.156"	69.088 mm 2.72"	10/100, 2 lb
	89031-284	1203-PFB BL VWR	VWR® Elongated tip, glass-filled pointed polypropylene handle	25.4 x 4.064 mm 1 x 0.16"	82.804 mm 3.26"	2/10/50, 2 lb
	89031-282	1223-PFB BL VWR	VWR® Reticulated foam cone shaped tip, glass-filled polypropylene handle	17.47 x 3.96 mm 0.688 x 0.156"	82.804 mm 3.26"	2/10/50, 2 lb
	89133-786	1273 BK	Soft flexible paddle tip, polypropylene handle	15.24 x 3.962 mm 0.6 x 0.156"	68.58 mm 2.7"	2/500, 1 lb
	89176-450	1273 BK VWR	VWR® Soft flexible paddle tip, black polypropylene handle			5/500, 2 lb
	89133-788	1276 BK	Soft flexible paddle tip, polypropylene handle	18.542 x 6.35 mm 0.73 x 0.25"	75.413 mm 2.969"	2/500, 4 lb
	89031-294	1276-PFB BL VWR	VWR® Soft flexible paddle tip, black polypropylene handle			5/500, 8 lb
	89031-286	1605-PSF RECT VWR	VWR® Rectangular tip, polypropylene handle	25.4 x 13.487 mm 1 x 0.531"	130.175 mm 5.125"	2/10/50, 5 lb
	89133-782	1616-PF	Standard tip, polypropylene handle	22.225 x 6.35 mm 0.875 x 0.25"	150.8 mm 5.937"	20/50, 3 lb
	89133-780	1802-PF SOLID	Standard tip, polypropylene handle	15.088 x 4.775 mm 0.594 x 0.188"	60.325 mm 2.375"	20/50, 2 lb
	10805-138	1803-PFB	Mini tip, polypropylene handle	15.875 x 3.962 mm 0.625 x 0.156"	71.45 mm 2.813"	2/500, 1 lb (bulk)
	89133-778	1803-WF	Standard tip, wood handle	11.913 x 3.962 mm 0.469 x 0.156"	74.625 mm 2.938"	20/50, 1 lb
	10806-020	1804-PF	Standard tip, polypropylene handle	15.875 x 3.962 mm 0.625 x 0.156"	101.6 mm 4"	20/50, 1 lb
	10806-021	1804-PFB				4/250, 2 lb
	10806-032	1805-PF RECT	Large rectangular tip, polypropylene handle	22.606 X 17.45 mm 0.89 X 0.687"	131.369 mm 5.172"	20/50, 6 lb
	10806-023	1805-PFB RECT				4/250, 6 lb (bulk)
	89133-776	1805-PF RND	Large round tip, polypropylene handle	28.575 x 23.012 mm 1.125 x 0.906"	130.175 mm 5.125"	20/50, 6 lb

PURSWAB® FOAM-TIPPED APPLICATORS, CONTINUED

	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	10806-024	1806-PF	Rectangular tip, polypropylene handle	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"	20/50, 2 lb
	10806-025	1806-PFB				4/250, 6 lb
	27500-100	1806-WF	Elongated tip, wood handle	17.45 x 4.775 mm 0.687 x 0.188"	149.225 mm 5.875"	20/50, 2 lb
	89176-730	1806-WFB				4/250, 2 lb
	89031-274	1806-WF VWR				VWR® Elongated tip, wood handle
	89031-276	1876-PF VWR	VWR® Flexible paddle tip, polypropylene handle	15.088 x 4.775 mm 0.594 x 0.188"	153.975 mm 6.062"	2/10/50, 3 lb

PURSWAB® FOAM OVER COTTON-TIPPED APPLICATORS

These foam over cotton tips are made from shred-resistant urethane foam over USP pharmaceutical-grade cotton for the best absorption and are great when cleaning with solvents. It's like using 2 swabs in 1!

	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	10806-026	1806-PCF	Foam-over-cotton tip, polypropylene handle	23.8 x 9.525 mm 0.937 x 0.375"	155.575 mm 6.125"	20/50, 2 lb
	10806-027	1806-PCFB				4/250, 2 lb
	10806-028	1806-PCFL	Large foam-over-cotton tip, polypropylene handle	25.4 x 9.525 mm 1 x 0.375"	155.575 mm 6.125"	20/50, 2 lb
	10806-152	1806-PCSF	Small foam-over-cotton tip, polypropylene handle	20.625 x 6.35 mm 0.812 x 0.25"	154.78 mm 6.094"	20/50, 2 lb
	82030-590	1806-WCF	Standard foam-over-cotton tip, wood handle	23.8 x 9.525 mm 0.937 x 0.375"	153.975 mm 6.062"	20/50, 2 lb
	89031-280	1806-WCF VWR				VWR® Standard foam-over-cotton tip, wood handle
	89031-278	1806-WCSF VWR	VWR® Small foam-over-cotton tip, wood handle	20.625 x 6.35 mm 0.812 x 0.25"	152.4 mm 6"	2/10/50, 2 lb
	82030-594	1806-WCSFB	Small foam-over-cotton tip, wood handle	20.625 x 6.35 mm 0.812 x 0.25"	152.4 mm 6"	4/250, 2 lb





PURSWAB® ANTI-STATIC APPLICATORS

Where minimizing electrostatic buildup is the goal, choose one of our anti-static-tipped applicators.

	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	10806-154	1806-PF AS	Standard anti-static foam tip, conductive polypropylene handle	15.875 x 3.962 mm 0.625 x 0.156"	150.8 mm 5.937"	20/50, 3 lb




PURSWAB® ESD APPLICATORS

These applicators are designed for working in and around static-sensitive components and parts. A range of shapes, sizes, and materials let you find a swab made specifically for each of your applications.

	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	89087-858	1273-PF ESD	Paddle anti-static foam tip, static-dissipative polypropylene handle	14.224 x 3.962 mm 0.56 x 0.156"	68.58 mm 2.7"	20/50, 2 lb
	89108-136	1283-PF ESD	Precision pointed foam tip, static-dissipative polypropylene handle	19.989 x 3.505 mm 0.787 x 0.138"	69.215 mm 2.725"	20/50, 1 lb
	89087-862	1621-PF ESD	Precision pointed foam tip, static-dissipative polypropylene handle	10.922 x 3.556 mm 0.43 x 0.14"	82.804 mm 3.26"	20/50, 2 lb
	89087-852	3605 ESD	Microfiber tip, static-dissipative polypropylene handle	6.129 x 5.996 mm 0.635 x 0.236"	152.4 mm 6"	10/100, 3 lb

PURSWAB® SPECIALTY + CHEMICAL RESISTANT FOAM-TIPPED APPLICATORS

This line of specialty swabs includes super-absorbent and chemical-resistant foam-tipped applicators. Chemical-resistant foam is made of a different polymer and will stand up to most solvents.

	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	10806-144	1806-WCF CR	Standard chemical resistant foam-over-cotton tip, wood handle	20.625 x 9.525 mm 0.812 x 0.375"	153.975 mm 6.062"	20/50, 2 lb
	10806-142	1806-WF CR	Standard chemical resistant foam tip, wood handle	17.45 x 4.775 mm 0.687 x 0.188"	149.225 mm 5.875"	20/50, 2 lb
	10806-146	1876-PF CR	Flexible chemical resistant paddle foam tip, polypropylene handle	15.088 x 4.775 mm 0.594 x 0.188"	153.975 mm 6.062"	20/50, 3 lb



PURITAN COTTON-TIPPED APPLICATORS

Ideal for general purposes when you are less concerned about static or particulate release.

	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	10805-140	803-PC	Standard tip, polystyrene handle	17.45 x 5.156 mm 0.687 x 0.203"	77.775 mm 3.062"	10/10/100, 9 lb
	10805-142	803-WC	Standard tip, wood handle	15.875 x 4.775 mm 0.625 x 0.188"	77.013 mm 3.032"	10/10/100, 10 lb
	10805-158	803-WCL	Large tip, wood handle	17.45 x 5.944 mm 0.687 x 0.234"	76.987 mm 3.031"	10/5/100, 7 lb
	10805-148	806-PC	Standard tip, polystyrene handle	17.45 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"	10/10/100, 17 lb
	10806-000	806-WC	Standard tip, wood handle	15.875 x 4.775 mm 0.625 x 0.188"	150.8 mm 5.937"	10/10/100, 16 lb
	89031-270	806-WC VWR	VWR® Standard tip, wood handle			10/10/100, 14 lb
	10806-001	806-WCL	Large tip, wood handle	17.45 x 5.944 mm 0.687 x 0.234"	150.8 mm 5.937"	10/5/100, 10 lb
	89031-272	806-WCL VWR	VWR® Large tip, wood handle			10/5/100, 10 lb
	89133-810	821-WC DBL	Double-tipped (one standard, one tapered), wood handle	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"	10/5/100, 9 lb
	10805-130	826-WC	Tapered mini tip, wood handle	12.7 x 3.429 mm 0.5 x 0.135"	149.225 mm 5.875"	10/10/100, 12 lb
	89176-684	836-WC NO BINDER	Standard tip (no binder), wood handle	15.875 x 4.77 mm 0.625 x 0.188"	77.01 mm 3.032"	10/10/100, 10 lb
	10806-008	860-WC	Standard tip, lint-sensitive, wood handle	15.875 x 4.775 mm 0.625 x 0.188"	150.8 mm 5.937"	10/10/100, 16 lb
	10806-004	861-WC DBL	Double standard tip, wood handle	15.875 x 4.775 mm 0.625 x 0.188"	153.975 mm 6.062"	10/5/100, 10 lb
	10806-005	867-WC NO GLUE	Standard tip (glue-free), wood handle	15.875 x 4.775 mm 0.625 x 0.188"	150.8 mm 5.937"	10/10/100, 16 lb
	82030-716	868-PCS	Small tip, polystyrene handle	17.45 x 4.369 mm 0.687 x 0.172"	152.4 mm 6"	10/10/100, 16 lb
	89237-580	868-WCS	Small tip, wood handle	11.1 x 3.581 mm 0.437 x 0.141"	149.225 mm 5.875"	10/10/100, 12 lb
	10806-010	869-WCS	Mini tip, wood handle	12.7 x 3.175 mm 0.5 x 0.125"	149.225 mm 5.875"	10/10/100, 13 lb

PURSWAB® LINT-SENSITIVE COTTON-TIPPED APPLICATORS

An excellent choice when you are looking for a swab for cleaning sensitive components. These swabs are both low-lint and residue-free.


	VWR Cat. No.	Puritan Item No.	Product Description	Tip length x diameter	Full Length	Packaging pack, weight
	82030-740	870-PC DBL	Double sharp pointed tip, fine paper handle	8 x 2.21 mm 0.315 x 0.087"	76.2 mm 3"	100/25, 2 lb
	82029-064	871-PC DBL	Double pointed tip, fine paper handle	10.998 x 2.997 mm 0.433 x 0.118"	78.994 mm 3.11"	100/25, 2 lb
	82030-744	872-PC DBL	Double cone-shaped tip, fine paper handle	10.998 x 4.7 mm 0.433 x 0.185"	81 mm 3.189"	50/25, 2 lb
	89031-288	872-PC DBL VWR	VWR® Double cone-shaped tip, fine paper handle			50/25, 2 lb
	89133-812	876-WC	Standard tip, lint-sensitive, wood handle	15.875 x 4.775 mm 0.625 x 0.188"	150.8 mm 5.937"	10/10/100, 15 lb
	82030-746	890-PC DBL	Double mini tip, fine paper handle	8.001 x 2.007 mm 0.315 x 0.079"	75.489 mm 2.972"	100/25, 2 lb
	82030-750	892-PC DBL	Double small tip, fine paper handle, anti-static packaging	13.995 x 3.2 mm 0.551 x 0.126"	75.489 mm 2.972"	100/25, 2 lb
	89031-292	892-PC DBL VWR	VWR® Double small tip, fine paper handle	13.995 x 3.2 mm 0.551 x 0.126"	75.489 mm 2.972"	100/25, 2 lb
	89031-290	893-PC DBL VWR	VWR® Double tip, paper handle, anti-static packaging	13.995 x 3.2 mm 0.551 x 0.126"	75.489 mm 2.972"	2/50/25, 4 lb
	82030-752	896-PC	Small cylindrical tip, paper handle, anti-static packaging	12 x 4 mm 0.47 x 0.157"	146 mm 5.748"	20/50, 3 lb
	82030-754	896-WC	Standard tip, wood handle	15 x 4.7 mm 0.59 x 0.185"	150 mm 5.905"	10/100, 2 lb

Did you know?

Only 25% of the North American hardwood processed becomes finished goods. That means only the “finest” part of the tree becomes tongue depressors, applicator sticks and other wood products. The rest becomes environmentally correct fuel for our operations.

INDUSTRIAL WOODEN STIRRERS


These are excellent lower-price alternatives to our medical-quality wood products. Ideal in manufacturing environments as they are consistently made to our tough standards.

	VWR Cat. No.	Puritan Item No.	Product Description	Full Length	Packaging pack, weight
	10805-001	704 BRIGHTWOOD	Flat blade, economy-grade, Brightwood® brand	152.4 mm 6"	10/500, 36 lb

WOODEN ORANGE STICKS

Made of North American hardwood, these were and often still are called orange wood sticks, as they were originally made of orange wood. Very handy, multi-purpose orange sticks are available with different end shapes.



	VWR Cat. No.	Puritan Item No.	Product Description	Dimensions length x diameter	Packaging pack, weight
	87004-104	2910	Double bevel, both ends	177.8 x 3.962 mm 7 x 0.156"	10/100, 5 lb

Puritan® – a tradition of excellence

Order with confidence, knowing you will always get consistent quality, ready availability, and outstanding service. You can count on Puritan.

From our inception in 1919, our wood products operation has grown to become a state-of-the-art medical device manufacturer serving markets worldwide with Puritan's US-made products.

Our success is evident in the success of our customers. Our guarantee is to produce the highest-quality single-use products possible. Our mission is to provide our customers with unparalleled service. Our promise is to deliver – the right product, on time.

Puritan Medical Products Co.
31 School Street
P.O. Box 149
Guilford, Maine 04443-0149 USA

800.321.2313 (US + Canada)
207.876.3311



Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2023 Avantor, Inc. All rights reserved.