



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Revision Date 01/26/2015

Version 1.2

## SECTION 1. Identification

### Product identifier

Product number HX0603  
Product name Hydrochloric Acid GR ACS

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses Reagent for analysis

### Details of the supplier of the safety data sheet

Company EMD Millipore Corporation | 290 Concord Road, Billerica, MA 01821,  
United States of America | General Inquiries: +1-978-715-4321 |  
Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)

Emergency telephone 800-424-9300 CHEMTREC (USA)  
+1-703-527-3887 CHEMTREC (International)  
24 Hours/day; 7 Days/week

## SECTION 2. Hazards identification

### GHS Classification

Corrosive to Metals, Category 1, H290  
Skin corrosion, Category 1B, H314  
Serious eye damage, Category 1, H318  
Specific target organ systemic toxicity - single exposure, Category 3, Respiratory system, H335  
For the full text of the H-Statements mentioned in this Section, see Section 16.

### GHS-Labeling

#### Hazard pictograms



#### Signal Word

Danger

#### Hazard Statements

H290 May be corrosive to metals.  
H314 Causes severe skin burns and eye damage.  
H335 May cause respiratory irritation.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number HX0603  
Product name Hydrochloric Acid GR ACS

Version 1.2

## *Precautionary Statements*

P234 Keep only in original container.  
P261 Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray.  
P264 Wash skin thoroughly after handling.  
P271 Use only outdoors or in a well-ventilated area.  
P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.  
P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
P303 + P361 + P353 IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower.  
P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P310 Immediately call a POISON CENTER or doctor/ physician.  
P321 Specific treatment (see supplemental first aid instructions on this label).  
P363 Wash contaminated clothing before reuse.  
P390 Absorb spillage to prevent material damage.  
P403 + P233 Store in a well-ventilated place. Keep container tightly closed.  
P405 Store locked up.  
P406 Store in corrosive resistant stainless steel container with a resistant inliner.  
P501 Dispose of contents/ container to an approved waste disposal plant.

## **Other hazards**

None known.

---

## **SECTION 3. Composition/information on ingredients**

Chemical nature Aqueous solution

### **Hazardous ingredients**

*Chemical Name (Concentration)*

CAS-No.

*hydrochloric acid (>= 30 % - < 50 % )*

7647-01-0

Exact percentages are being withheld as a trade secret.

---

## **SECTION 4. First aid measures**

### **Description of first-aid measures**

*General advice*

First aider needs to protect himself.

*Inhalation*

After inhalation: fresh air. Call in physician.

*Skin contact*

In case of skin contact: Take off immediately all contaminated clothing. Rinse skin with water/ shower. Call a physician immediately.

*Eye contact*

After eye contact: rinse out with plenty of water. Immediately call in ophthalmologist.

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number HX0603  
Product name Hydrochloric Acid GR ACS

Version 1.2

### *Ingestion*

After swallowing: make victim drink water (two glasses at most), avoid vomiting (risk of perforation!). Call a physician immediately. Do not attempt to neutralize.

Never give anything by mouth to an unconscious person.

### **Most important symptoms and effects, both acute and delayed**

Irritation and corrosion, Cough, Shortness of breath, cardiovascular disorders, Risk of blindness!

### **Indication of any immediate medical attention and special treatment needed**

No information available.

---

## **SECTION 5. Fire-fighting measures**

### **Extinguishing media**

#### *Suitable extinguishing media*

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### **Special hazards arising from the substance or mixture**

Not combustible.

Ambient fire may liberate hazardous vapors.

Fire may cause evolution of:

Hydrogen chloride gas

### **Advice for firefighters**

#### *Special protective equipment for fire-fighters*

Stay in danger area only with self-contained breathing apparatus. Prevent skin contact by keeping a safe distance or by wearing suitable protective clothing.

#### *Further information*

Suppress (knock down) gases/vapors/mists with a water spray jet. Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## **SECTION 6. Accidental release measures**

### **Personal precautions, protective equipment and emergency procedures**

Advice for non-emergency personnel: Do not breathe vapors, aerosols. Avoid substance contact. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### **Environmental precautions**

Do not empty into drains.

### **Methods and materials for containment and cleaning up**

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up with liquid-absorbent and neutralizing material (e.g. Chemisorb® H<sup>+</sup>, Art. No. 101595).

Dispose of properly. Clean up affected area.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

HX0603

Version 1.2

Product name

Hydrochloric Acid GR ACS

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

### Conditions for safe storage, including any incompatibilities

*Requirements for storage areas and containers*

No metal containers.

Tightly closed.

Store at room temperature.

---

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

#### *Ingredients*

Basis	Value	Threshold limits	Remarks
<i>hydrochloric acid 7647-01-0</i>			
ACGIH	Ceiling Limit Value:	2 ppm	
NIOSH/GUIDE	Ceiling Limit Value and Time Period (if specified):	5 ppm 7 mg/m <sup>3</sup>	
OSHA_TRANS	Ceiling Limit Value:	5 ppm 7 mg/m <sup>3</sup>	
Z1A	Ceiling Limit Value:	5 ppm 7 mg/m <sup>3</sup>	

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

#### *Hygiene measures*

Immediately change contaminated clothing. Apply skin- protective barrier cream. Wash hands and face after working with substance.

#### *Eye/face protection*

Tightly fitting safety goggles

#### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

#### *Other protective equipment:*

Acid-resistant protective clothing.

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

HX0603

Version 1.2

Product name

Hydrochloric Acid GR ACS

---

### *Respiratory protection*

required when vapors/aerosols are generated.

Use a properly fitted, air-purifying or air-fed respirator complying with an approved standard if a risk assessment indicates this is necessary. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator.

---

## SECTION 9. Physical and chemical properties

Physical state	liquid
Color	colorless
Odor	stinging
Odor Threshold	0.8 - 5 ppm Gaseous hydrogen chloride (HCl).
pH	< 1 at 68 °F (20 °C)
Solidification point	-30 °C
Boiling point	No information available.
Flash point	Not applicable
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	Not applicable
Upper explosion limit	Not applicable
Vapor pressure	190 hPa at 68 °F (20 °C)
Relative vapor density	No information available.
Density	ca. 1.19 g/cm <sup>3</sup> at 68 °F (20 °C)
Relative density	No information available.
Water solubility	at 68 °F (20 °C) soluble
Partition coefficient: n-octanol/water	Not applicable

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

HX0603

Version 1.2

Product name

Hydrochloric Acid GR ACS

---

Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	2.3 mPa.s at 59 °F (15 °C)
Explosive properties	Not classified as explosive.
Oxidizing properties	none
Ignition temperature	Not applicable
Corrosion	May be corrosive to metals.

---

## SECTION 10. Stability and reactivity

### Reactivity

Corrosive in contact with metals

### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

### Possibility of hazardous reactions

Exothermic reaction with:

Amines, potassium permanganate, salts of oxyhalogenic acids, semimetallic oxides, semimetallic hydrogen compounds, Aldehydes, vinylmethyl ether

Risk of ignition or formation of inflammable gases or vapors with:

carbides, lithium silicide, Fluorine

Generates dangerous gases or fumes in contact with:

Aluminum, hydrides, formaldehyde, Metals, strong alkalis, Sulfides

Risk of explosion with:

Alkali metals, conc. sulfuric acid

### Conditions to avoid

Heating.

### Incompatible materials

Metals, metal alloys

Gives off hydrogen by reaction with metals.

### Hazardous decomposition products

in the event of fire: See section 5.

---

## SECTION 11. Toxicological information

### Information on toxicological effects

*Likely route of exposure*

Inhalation, Eye contact, Skin contact

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

HX0603

Version 1.2

Product name

Hydrochloric Acid GR ACS

---

### *Target Organs*

Eyes

Skin

Respiratory system

Cornea

### *Acute oral toxicity*

Symptoms: If ingested, severe burns of the mouth and throat, as well as a danger of perforation of the esophagus and the stomach.

Acute toxicity estimate: 1,892 mg/kg

Calculation method

### *Acute inhalation toxicity*

Symptoms: mucosal irritations, Cough, Shortness of breath, Possible damages:, damage of respiratory tract

Acute toxicity estimate: 6.41 mg/l; 4 h

Calculation method

### *Skin irritation*

Mixture causes burns.

### *Eye irritation*

Mixture causes serious eye damage. Risk of blindness!

### *Specific target organ systemic toxicity - single exposure*

Target Organs: Respiratory system

Mixture may cause respiratory irritation.

### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

## **Carcinogenicity**

IARC

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

OSHA

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

NTP

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

ACGIH

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

## **Further information**

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

HX0603

Version 1.2

Product name

Hydrochloric Acid GR ACS

---

After uptake:

After a latency period:

cardiovascular disorders

Handle in accordance with good industrial hygiene and safety practice.

## Ingredients

*hydrochloric acid*

No information available.

---

## SECTION 12. Ecological information

### Ecotoxicity

No information available.

### Persistence and degradability

No information available.

### Bioaccumulative potential

*Partition coefficient: n-octanol/water*

Not applicable

### Mobility in soil

No information available.

### *Additional ecological information*

Forms corrosive mixtures with water even if diluted. Harmful effect due to pH shift.

Discharge into the environment must be avoided.

## Ingredients

*hydrochloric acid*

Substance does not meet the criteria for PBT or vPvB according to Regulation (EC) No 1907/2006, Annex XIII.

---

## SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.



SAFETY DATA SHEET  
according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number HX0603  
Product name Hydrochloric Acid GR ACS

Version 1.2

---

**SECTION 14. Transport information**

**Land transport (DOT)**

UN number UN 1789  
Proper shipping name HYDROCHLORIC ACID  
Class 8  
Packing group II  
Environmentally hazardous --

**Air transport (IATA)**

UN number UN 1789  
Proper shipping name HYDROCHLORIC ACID  
Class 8  
Packing group II  
Environmentally hazardous --  
Special precautions for user no

**Sea transport (IMDG)**

UN number UN 1789  
Proper shipping name HYDROCHLORIC ACID  
Class 8  
Packing group II  
Environmentally hazardous --  
Special precautions for user yes  
EmS F-A S-B

---

**SECTION 15. Regulatory information**

**United States of America**

**SARA 313**

The following components are subject to reporting levels established by SARA Title III, Section 313:

*Ingredients*

hydrochloric acid 7647-01-0 37 %

**SARA 302**

The following components are subject to reporting levels established by SARA Title III, Section 302:

*Ingredients*

hydrochloric acid 7647-01-0

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

HX0603

Version 1.2

Product name

Hydrochloric Acid GR ACS

---

## Clean Water Act

The following Hazardous Substances are listed under the U.S. CleanWater Act, Section 311, Table 116.4A:

### *Ingredients*

hydrochloric acid

The following Hazardous Chemicals are listed under the U.S. CleanWater Act, Section 311, Table 117.3:

### *Ingredients*

hydrochloric acid

## DEA List I

Not listed

## DEA List II

Listed

### *Ingredients*

hydrochloric acid

7647-01-0

## US State Regulations

### Massachusetts Right To Know

#### *Ingredients*

hydrochloric acid

### Pennsylvania Right To Know

#### *Ingredients*

hydrochloric acid

### New Jersey Right To Know

#### *Ingredients*

hydrochloric acid

### California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

## Notification status

TSCA: All components of the product are listed in the TSCA-inventory.

DSL: All components of this product are on the Canadian DSL.

KOREA: Not in compliance with the inventory

---

## SECTION 16. Other information

### Training advice

Provide adequate information, instruction and training for operators.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number

HX0603

Version 1.2

Product name

Hydrochloric Acid GR ACS

---

## Labeling

### *Hazard pictograms*



### *Signal Word*

Danger

### *Hazard Statements*

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

H335 May cause respiratory irritation.

### *Precautionary Statements*

#### Prevention

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

#### Response

P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308 + P310 IF exposed or concerned: immediately call a POISON CENTER or doctor/ physician.

### **Full text of H-Statements referred to under sections 2 and 3.**

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
H318	Causes serious eye damage.
H335	May cause respiratory irritation.

### **Key or legend to abbreviations and acronyms used in the safety data sheet**

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Revision Date 01/26/2015

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M" Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*