



# Product Specification

**Catalog Number:** BDH9278

**Catalog Description:** BDH SODIUM ACETATE ACS

**Grade:** ACS

**Date Current Specification Set:** NOV. 2012

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

Characteristics	Specifications
Calcium	<= 0.005 %
Chloride	<= 0.001 %
Heavy Metals (as Pb)	<= 0.0005 %
Insolubles	<= 0.005 %
Iron	<= 0.0005 %
Magnesium	<= 0.002 %
Permanganate Reducers Test	PASS
pH (5%, Water) @25C	7.5 - 9.2
Phosphate	<= 0.0005 %
Potassium	<= 0.005 %
Purity	99.0 - 101.0 %
Sulfate	<= 0.002 %

Signed on behalf of VWR:

Robin Crocker  
Quality Control Manager