

BULK DURABLE TRANSPORT

PROTECT PRODUCT QUALITY WITH BULK DURABLE INSULATED CONTAINERS

ThermoSafe® Durable Insulated Containers keep your temperature-sensitive products at the proper temperature and preserve product quality. Use these bulk insulated containers to optimize mixed loads of refrigerated, frozen, and controlled room temperature product. Our durable insulated containers are also ideal for temporary storage in the warehouse, on the pack out line, and for the storage refrigerants and dry ice.

Ranging in various sizes and styles, these durable insulated containers make storing, transporting, and unloading efficient and worry-free. High impact polyethylene can handle the rigors of transport: the tight fitting lid and container construction, maintains temperatures, and maximizes shipping and transporting efficiency.




TOP LOAD MODELS

...IMAGINE THE POSSIBILITIES

Constructed from FDA approved materials and ergonomically designed to promote operational efficiencies and to maximize cube space within a truck while maintaining the required temperature.

- Double-wall construction with specialty treated high impact polyethylene enhances foam adhesion and adds structural strength
- Specially formulated foam-in-place urethane insulation provides outstanding temperature control and long lasting durability
- Available in top load or front load style ranging from 4 to 74 cubic feet capacity, and can hold up to 1450 lb of pellets or cut dry ice
- FLX-SEAL® doors on upright models fold up and out of the way, allowing for complete product access
- Optional accessories include: dry ice compartments, shelving, drains, casters, airflow floors, retention straps



MODEL	DIMENSIONS L x W x H (in)		EMPTY WEIGHT (lb)	CAPACITY			
	INNER	OUTER		CUBIC FEET	DRY ICE		AVG. HOLDING TIME PELLETS*
					PELLETS (lb)	BLOCKS	
89496-032	22 x 12 x 26	26 x 16 x 33	47	4	200	4	12 days
15715-168	36½ x 21 x 24	43 x 27½ x 35¼	106	11	550	12	18 days
87003-916	36½ x 21 x 24	43 x 27½ x 39½	131	11	550	12	18 days
10818-628	37 x 41 x 20	43 x 48 x 31¾	165	19	990	24	21 days
10818-630	37 x 41 x 29	43 x 48 x 40½	188	27	1,400	36	25 days
10818-632	37 x 41 x 31	43 x 48 x 41½	193	30	1,450	36	25 days
10818-634	37 x 41 x 31	43 x 48 x 41½	243	30	1,450	36	25 days
10818-636	54 x 41 x 68½	61¼ x 48½ x 82	550	74	N/A	N/A	N/A



FRONT LOAD MODELS

ThermoSafe chests and containers are strong, durable, and specifically manufactured to be used and reused. Replacement latches, hinges, handles, and other handling and replacement accessories are available at a low cost. Contact your VWR Sales Representative for more information.



1.800.932.5000
vwr.com



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2016 VWR International, LLC. All rights reserved.