

## VWR<sup>®</sup> REPEATING DISPENSER & SYRINGE TIPS

Light-weight,  
ergonomic design

Meets DIN EN ISO  
8655 requirements

Allows for nine tip  
sizes, ten volume  
adjustments,  
and ninety  
program steps



## VWR® Repeating Dispenser

Manufactured under DIN EN ISO 8655 requirements using high-quality materials, ensuring a device that is built to last. Volume selection, loading, and dispensing of liquids can be performed with ease using only one hand. The smooth-moving dosing button and light-weight design facilitates long dosing series without causing strain or discomfort to the user.

- Each device is tested for functionality and precision
- Maintenance-free; does not require batteries or a power supply
- Resistant to common laboratory chemicals
- Light-weight, ergonomic design
- Volume adjustments and dosing volumes are easily visible in panel at top of pipettor
- Designed to precisely align with VWR® Classic Syringe Tips
- Compatible with Eppendorf's Combitips®, Combitips Plus®, and Combitips Advanced®



Description	Cat. No.	Unit
VWR Repeating Dispenser	75836-686	Each
<b>Accessories</b>		
Adapter for 25 mL Pro Style Dispenser Tips	75836-688	Each
Adapter for 50 mL Pro Style Dispenser Tips	75836-690	Each

## VWR® Classic Syringe Tips

VWR Classic Syringe Tips are ideal for viscous, infectious, and volatile liquids, and are compatible with a variety of pipettors. Sterile tips are individually wrapped, and certified to be pyrogen-, ATP-, RNase-, and DNase-free.

Volume, mL	Cat. No.	Pack of
<b>VWR Classic Syringe Tips, Non Sterile</b>		
0.05	76163-558	100
0.5	76163-560	100
1.25	76163-592	100
2.5	76163-594	100
5.0	76163-596	100
12.5	76163-598	100
25.0	75838-208	25
50.0	75838-210	25
<b>VWR Classic Syringe Tips, Sterile</b>		
0.05	76163-600	100
0.5	76163-602	100
1.25	76163-604	100
2.5	76163-606	100
5.0	76163-608	100
12.5	76163-610	100
25.0	75838-226	25
50.0	75838-228	25



Suitable for use with the following pipettors: Ripette®, Ripette® pro, Ripette® genX, eLine™-Dispenser, Minilab 100, Minilab 101, Minilab 201, Distriman®, StepMate, EasyStep, EasyStep electronic, RepeatOne®, and Handrop.