

Product Name:	VWR-HCSSP-49FA
Description:	<p>VWR Catalog Number: 76514-218 VWR Plus Series Stainless Steel Freezer 49 Cu. Ft. capacity Two swing solid doors, self-closing, not reversible 6 Ea. Adjustable Shelves Digital display microprocessor temperature controller Adjustable operating temperature: -15°C to -25°C Audible and visual Hi/Lo temperature alarm 2 probes; 1 in air and 1 in sample bottle Probe attached to evaporator for defrost termination sensing Sensor error alarm Continuous product (simulated) temperature display Remote alarm contacts Forced draft air circulation Keyed door lock Magnetic door gaskets for positive seal Integrated door handle LED interior lights are safety shielded and switch controlled Probe access port (3/4") Stainless steel exterior High density urethane foam cabinet insulation White-powder coated aluminum interior with 304 stainless steel floor Easily cleaned smooth interior Hydrocarbon natural refrigerant (R290) Two year parts and labor warranty, plus an additional three year compressor parts warranty Auto defrost 4 swivel casters (5 3/4") 115V/60 Hz ETL listed Exterior dimensions: 54" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters 115V, 60 Hz, 9 Amps, 3/4 HP ETL listed Shipping weight: 495 lbs.</p>