

## Anti-DYKDDDDK conjugated to Dylight® 650

Product Number D10-1718  
Amount 100 µg or 1 mg total protein  
Store at 2-8°C

### Form/Shipping & Storage

Supplied as a liquid. Upon receipt, store at 2-8 °C. Fluorescent dyes are sensitive to light. Please store in the dark and avoid exposure to sunlight.

### Handling

We recommend that the investigator determine the appropriate working concentration for their specific application.

### Buffer

0.1M sodium phosphate, 0.15M NaCl, pH7.4, 1% BSA, 0.02% NaN<sub>3</sub>

### Stability


Product should be stored at 2-8°C in the dark and be used within 1 year. If further dilution of the conjugate is required, use diluted material within one month.

### Note

For research use only, not for diagnostic or therapeutic use.

### Spectral Characteristic

Visible absorption maxima 650 nm  
Emission maximum 672 nm



### References

Puertollano R, Alonso MA. MAL, an integral element of the apical sorting machinery, is an itinerant protein that cycles between the trans-Golgi network and the plasma membrane. *Mol Biol Cell*. 1999 Oct;10(10):3435-3447.

Roy SF, Laporte SA, Escher E, Leduc R, Guillemette G. Epitope tagging and immunoreactivity of the human angiotensin II type 1 receptor. *Can J Physiol Pharmacol*. 1997 Jun;75(6):690-695.



**COLUMBIA BIOSCIENCES**