

L-Shape Cell Spreader

Description: Purple, Sterile, non-pyrogenic


Purpose: Disposable L-shaped cell spreaders are ideal for spreading cell or bacterial cultures evenly

Materials

Cell Spreader: PP (Polypropylene)

Color: Purple

Dimensions (Unit: mm)



Features

- Sterilization: sterilized by gamma irradiation SAL10⁶(ISO11137)
- DNase/RNase free, non-pyrogenic
- Individually wrapped and bulk packaging
- Disposable L-shaped cell spreaders are ideal for spreading cell or bacterial cultures evenly
- The L-shaped allows application on the entire dish or plate surface
- Smooth spreading surface with slight upward return reduces the risk of damaging agar sample
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under EN ISO13485:2003 and ISO9001:2008 quality management system

L-Shape Cell Spreader

VWR EU EAN. No.	VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
	76208-438	Purple, Sterile, Non-pyrogenic, Individually wrapped	Yes	1	100	100
	76207-748	Purple, Sterile, Non-pyrogenic, Bulk Package	Yes		10	500