



**Item Number** C1977  
**Item** Calcium Chloride, Anhydrous, Pellets, FCC  
**CAS Number** 10043-52-4  
**Molecular Formula** CaCl<sub>2</sub>  
**Molecular Weight** 110.98

Assay (CaCl<sub>2</sub>).....93.0-100.5%

MAXIMUM LIMITS

Acid-Insoluble Matter..... 0.02%  
Arsenic (as As)..... 3 mg/kg  
Fluoride..... 0.004%  
Lead..... 5 mg/kg  
Magnesium and Alkali Salts..... 5.0%

***Spectrum Chemicals & Laboratory Products***

---

**Corporate Headquarters:**      **East Coast Facility:**  
14422 S. San Pedro St.      755 Jersey Ave.  
Gardena, CA 90248      New Brunswick, NJ 08901  
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843