

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Alliance Microporous Shoe covers, Universal shoe size, white color
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-03-21

VWR Catalogue Number	Description
76491-696	Alliance microporous shoecovers
76491-702	Alliance microporous shoecovers, anti-skid

Quality System Compliance

Products are manufactured in accordance with the manufacture's quality system. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	White alliance breathable microporous polyethylene coated polypropylene made in USA.
Sterilization	Product is non-sterile.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	This product is non-toxic.

Non-Pyrogenic Statement	Product contain polypropylene / polyethylene material.
Latex Statement	The product is Latex free.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Diethylhexylphthalate (DEHP) is not included in the manufacture or the formulation of the raw material and are not expected to be present.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store at ambient temperature.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.