



Material 76432-700
Material description Isopropyl alcohol 70%

CAS Number 67-63-0
Molecular formula CH3CHOHCH3
Molecular mass 60.10

Additional information

Characteristics

Specifications

Assay (V/V)	69.0 --> 71.0 %
Appearance	Clear liquid
Specific gravity	0.873 --> 0.883
Prepared to formulation on file	Confirmed

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada

Additional information

The weights and/or volumes used to prepare this product are NIST traceable. All balances used in the preparation of the product are calibrated weekly against NIST traceable weights. The balances are maintained and serviced on a regular basis by an outside certified company. All volumetric glassware used is NIST traceable and certified as meeting Class A specifications. Unless otherwise indicated above, all chemicals used in the preparation of this product are Reagent ACS grade.