

Catalog Number: BDH9244

Catalog Description: BDH MAGNESIUM CHLORIDE ACS

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics

Specifications

Ammonium	<= 0.002 %
Barium	<= 0.005 %
Calcium	<= 0.01 %
Heavy Metals (as Pb)	<= 0.0005 %
Insolubles	<= 0.005 %
Iron	<= 0.0005 %
Manganese	<= 0.0005 %
Nitrate	<= 0.001 %
Phosphate	<= 0.0005 %
Potassium	<= 0.005 %
Purity	99.0 - 102.0 %
Sodium	<= 0.005 %
Strontium	<= 0.005 %
Sulfate	<= 0.002 %

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager