

# ROTATOR GENIE™

Rugged and Versatile Lab Rotator!



- © Holds a wide variety of tubes, flasks, microplates, bags, and odd-shaped vessels at any angle.
- © Available in two INTERCHANGEABLE platform configurations.
- © MAGNETIC platform features unique spring clip plate system. Standard platform turns 90 degrees.
- © Designed for cold room and incubator use. Will easily fit in most incubators.
- © Uniquely designed for a wide selection of motions adjustable from gentle rolling to forceful end over end mixing.

## Specifications

Supplied	Standard Platform with Retaining Bands or Magnetic Platform
Speed	5-35 RPM
Overall Dimensions (W x D x H)	400 x 175 x 400mm (15.8 x 7.0 x 15.8 in) Depth with Platform Horizontal - 290mm (11.5 in)
Platform Dimensions	
Standard Platform	260 x 260mm (10.25 x 10.25 in)
Magnetic Platform	290 x 315mm (11.5 x 12.5 in)
Weight	5.2 kgs (11.5 lbs)
Environmental	Installation Category 2, 0° - 38°C (32° - 100°F), 95% max. relative humidity



Fig. No.	VWR Cat. No.	Description	Fig. No.	VWR Cat. No.	Description
--	<b>89202-302</b>	Rotator Genie, Standard Platform, 120V	6	<b>15001-002(*)</b>	6 Clip Plates for 12 each 10-13mm tubes
--	<b>89202-306</b>	Rotator Genie, Magnetic Platform, 120V	7	<b>15001-004(*)</b>	6 Clip Plates for 6 each 15-17mm tubes
1	<b>89202-310</b>	Standard Platform with Bands	8	<b>15001-006(*)</b>	6 Clip Plates for 3 each 28-30mm tubes
2	<b>89202-312</b>	Magnetic Platform	9	<b>15001-014(*)</b>	Bag Mounting Strips (pack of 16)
3	<b>15001-008(*)</b>	1 Clip Plate for 12 each 10-13mm tubes	10	<b>14003-576(*)</b>	2 Clip Plates for Hybridization Tubes
4	<b>15001-010(*)</b>	1 Clip Plate for 6 each 15-17mm tubes	11	<b>14222-140(*)</b>	2 Universal Clip Plates
5	<b>15001-012(*)</b>	1 Clip Plate for 3 each 28-30mm tubes	12	<b>14222-144(*)</b>	1 Magnetic Covered 100 Microtube Box

(\*) - Requires magnetic platform, (Cat. No. 89202-312), which is supplied with all 89202-306 models or can be purchased separately.



# ROTO-SHAKE GENIE

ROCK, ROTATE OR 3-D ORBIT ANY CONTAINER!



U.S. Pat. 5,918,979

- Platform ROCKS or ROTATES
- Magnetic Platform and clip plate system
- Optional Stackable Tray doubles ROCKING capacity
- May be used in cold rooms and incubators

### Specifications

Supplied	1 Clip Plate each for 10-13 mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray
Rotating	1-35 RPM
Rocking	2-70 cycles per minute
Platform Capacity	4.5 Kg (10 lb)
Weight	11.8 Kg (26lb)
Dimensions (D x W x H)	252 x 508 x 267mm
(Overall)	9.9 x 20 x 10.5in
Depth with Rocker Tray	305mm (12in)
Environmental	Installation Category 2, 0° - 38°C (32° - 100°F), 95% max. relative humidity
Warranty	2 Years



Fig.	VWR Cat. No.	Description	Fig.	VWR Cat. No.	Description
--	<b>58815-176</b>	Roto-Shake Genie, 120V	8	<b>15001-016</b>	Bag Expansion Kit
1	<b>15001-008</b>	Clip Plate for 12 each of 10-13mm tubes, 1/pk	9	<b>14003-576</b>	Clip Plates for Hybridization Tubes, 2/pk
2	<b>15001-010</b>	Clip Plate for 6 each of 15-17mm tubes, 1/pk	10	<b>89401-428</b>	Stackable Tray, 1/pk
3	<b>15001-012</b>	Clip Plate for 3 each of 28-30mm tubes, 1/pk	11	<b>14222-140</b>	Universal Clip Plates, 2/pk
4	<b>15001-002</b>	Clip Plates for 12 each of 10-13mm tubes, 6/pk	12	<b>14222-144</b>	Magnetic Attachable Covered 100 Microtube Box, 1/pk
5	<b>15001-004</b>	Clip Plates for 6 each of 15-17mm tubes, 6/pk	13	<b>14003-580</b>	3-D Orbital Shaker Attachment
6	<b>15001-006</b>	Clip Plates for 3 each of 28-30mm tubes, 6/pk	14	<b>10753-556</b>	Horizontal Orbital Shaker Attachment
7	<b>15001-014</b>	Bag Mounting Strips, 16/pk			

# ENVIRO-GENIE®



Accessories Supplied with Unit

The multi-function Enviro-Genie® is designed to chill, heat, rotate, rock and stir all in one compact unit.

- ☉ Magnetic stainless steel platform may be easily set to rock or rotate almost any container
- ☉ Containers can be easily mounted on the magnetic platform with magnetic clips or strips
- ☉ Optional 3-D orbital attachment transforms the unit into a temperature-controlled orbital platform shaker
- ☉ RS232 interface for data download and/or remote operation

## Specifications

Supplied	1 Clip Plate each for 10-13mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray, Shelf, 2 Integrated Magnetic Stirrers
Temperature	Heating: to 75°C (167°F), Cooling: to 4°C (39°F) Accuracy: ± 0.2°C, Uniformity: ± 0.5°C Resolution: 0.1°C, Stability: ± 0.1°C
Rotating	1-35 RPM
Rocking	2-70 cycles per minute
Stirring	1-2000 RPM (including Stir Reverse)
Platform Capacity	4.5 Kgs (10 lbs)
Chamber Volume	24 L (0.86 ft <sup>3</sup> )
Weight	36 Kgs (80 lbs)
Dimensions (D x W x H) (Overall)	495 x 622 x 368 mm
Dimensions (D x W x H) (Chamber)	19.5 x 24.5 x 14.5 in
Time, Hours	254 x 362 x 260 mm
Alarms	10 x 14.25 x 10.25 in Programmable to 99 Visible and Audible for Temperature and Time



Fig.	VWR Cat. No.	Description	Fig.	VWR Cat. No.	Description
--	<b>13500-846</b>	Enviro-Genie	7	<b>15001-014</b>	Bag Mounting Strips (pack of 16)
1	<b>15001-008</b>	1 Clip Plate for 12 each 10-13mm tubes	8	<b>15001-016</b>	Bag Expansion Kit
2	<b>15001-010</b>	1 Clip Plate for 6 each 15-17mm tubes	9	<b>14003-576</b>	2 Clip Plates for Hybridization Tubes
3	<b>15001-012</b>	1 Clip Plate for 3 each 28-30mm tubes	10	<b>14003-582</b>	1 Stackable Shelf
4	<b>15001-002</b>	6 Clip Plates for 12 each 10-13mm tubes	11	<b>14222-140</b>	2 Universal Clip Plates
5	<b>15001-004</b>	6 Clip Plates for 6 each 15-17mm tubes	12	<b>14222-144</b>	1 Magnetic Covered 100 Microtube Box
6	<b>15001-006</b>	6 Clip Plates for 3 each 28-30mm tubes	13	<b>14003-574</b>	3-D Orbital Shaker Attachment
			14	<b>490015-482</b>	Horizontal Orbital Shaker Attachment

# INCUBATOR-GENIE™

75.0°C  
28.0°C



Accessories Supplied with Unit



U.S. Pat. 5,918,979

- Ideal for hybridization, digestions, or incubation
- Magnetic stainless steel platform may be easily set to rock or rotate almost any container
- Containers can be easily mounted on the magnetic platform with magnetic clips or strips
- Chamber air is gently and continuously re-circulated ensuring temperature consistency and uniformity

### Specifications

Supplied	1 Clip Plate each for 10-13 mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray, Shelf
Temperature	28°C - 75°C (82°F - 167°F) Accuracy: ± 0.2°C, Uniformity: ± 0.5°C Resolution: 1.0°C, Stability: ± 0.1°C
Rotating	3-35 RPM
Rocking	6-70 cycles per minute
Platform Capacity	4.5 Kg (10 lb)
Chamber Volume	24 L (0.86 ft³)
Weight	19.5 Kg (43 lb)
Dimensions (D x W x H) (Overall)	400 x 560 x 370mm
Dimensions (D x W x H) (Chamber)	15.75 x 22 x 14.5in
Alarms	254 x 362 x 260mm
Warranty	10 x 14.25 x 10.25in Visible and Audible for Temperature 2 Years



Fig.	VWR Cat. No.	Description	Fig.	VWR Cat. No.	Description
--	<b>97060-138</b>	Incubator Genie, 120V	8	<b>15001-016</b>	Bag Expansion Kit
1	<b>15001-008</b>	Clip Plate for 12 each of 10-13mm tubes, 1/pk	9	<b>14003-576</b>	Clip Plates for Hybridization Tubes, 2/pk
2	<b>15001-010</b>	Clip Plate for 6 each of 15-17mm tubes, 1/pk	10	<b>14003-582</b>	Stackable Shelf, 1/pk
3	<b>15001-012</b>	Clip Plate for 3 each of 28-30mm tubes, 1/pk	11	<b>14222-140</b>	Universal Clip Plates, 2/pk
4	<b>15001-002</b>	Clip Plates for 12 each of 10-13mm tubes, 6/pk	12	<b>14222-144</b>	Magnetic Attachable Covered 100 Microtube Box, 1/pk
5	<b>15001-004</b>	Clip Plates for 6 each of 15-17mm tubes, 6/pk	13	<b>14003-574</b>	3-D Orbital Shaker Attachment
6	<b>15001-006</b>	Clip Plates for 3 each of 28-30mm tubes, 6/pk	14	<b>490015-482</b>	Horizontal Orbital Shaker Attachment
7	<b>15001-014</b>	Bag Mounting Strips, 16/pk			